

Technical data sheet

Rotavapor® R-250 Pro

Batches up to 30 liters can be quickly and safely distilled in the 50-liter evaporating flask of the R-250 Pro. The powerful heating bath allows a distillation rate up to 30 liters of acetone per hour. Choose your system out of a wide range of possibilities to get the best for your application.



Scope of delivery

All configurations are supplied ready to use.

| Components | R-250 Pro |
|--|-----------|
| Glass configuration according to order | 1 |
| Evaporating flask according to order | 1 |
| All needed tubing's (vacuum and cooling) | 1 |
| Seal removal tool | 1 |

Order code

Choose the configuration according to your needs:



Evaporating flask

D 50-liter evaporating flask

X no flask

Condensation unit (glass assembly)

51 R2 Double reflux

52 RB2 Double reflux Bullfrog

54 D2 Double descending





55 D3 Double descending with secondary condenser

Collection unit

1 Single receiving 20 liter

2 Dual receiving flasks 2 x 20 liters

Glassware



Assembly D2

Assembly D3

Assembly RB2

Assembly R2

| | | | |
|--|-------------------|---------------------|-------------|
| Low boiling points and/or foaming products | | High boiling points | |
| | Minimum emissions | Reflux reactions | |
| | | Reduced height | |
| H = 2300 mm | H = 2300 mm | H = 2100 mm | H = 2260 mm |

Technical data

Rotavapor® R-250 Pro

| | |
|---|--------------------------|
| Dimensions (W x D x H) (without glass) | 1410 x 830 x 1550 mm |
| Dimensions (W x D x H) (with glass) | 1450 x 850 x 2300 mm |
| Minimum clearance on all sides | 400 mm |
| Weight (without glass) | 160 kg |
| Weight (with glass) | 200 kg |
| Connection voltage | 400 ± 10 % VAC (3Ph,N,G) |
| Power consumption | 7500 W |
| Frequency | 50 / 60 Hz |
| IP Code | IP20 |
| Pollution degree | 2 |
| Overvoltage category | II |
| Pump outlet | max. 2 A |
| Rotation speed range | 5 – 120 rpm |
| Heating bath temperature range | 20 - 180 °C ± 2 °C |

| | |
|---|---|
| Adjustment accuracy | ± 1 °C |
| Regulation precision | At 60 °C: ± 1 °C At 95 °C: ± 2 °C At 180 °C: ± 3 °C |
| Heating medium | Water Polyethylene glycol 400 |
| Minimum flashpoint of the heating oil | 205 °C |
| Cooling water consumption | 200 - 400 L/h |
| Vacuum pump requirement | min 3 m ³ / h |
| Leakage of the complete system | <1 mbar/min |
| Approval (400 VAC Connection Voltage) | CE UL / CSA |
| Rotation controlling | Electronically |
| Rotation accuracy | ± 1 rpm at 5 rpm to ± 5 rpm at 120 rpm |
| Cooling restriction abs. without pulsation | max. 2.7 bar |
| Heating capacity | 6600 W |

Ambient conditions

For indoor use only.

| | |
|---------------------------------|--|
| Max. altitude above sea level | 2000 m |
| Ambient and storage temperature | 5 – 40 °C |
| Maximum relative humidity | 80% for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C |

Materials

| | |
|-------------------------|------------------------------------|
| Housing | Stainless steel 1.4301 (AISI 304) |
| Gear head | Aluminum cast (3.2373) |
| Painting | Powder coated with Epoxy (EPX) |
| Bath pan | Stainless steel 1.4404 (AISI 316L) |
| Heating element | Stainless steel 1.4404 (AISI 316L) |
| Glass | Borosilicate 3.3 |
| In contact with product | FDA approved materials |

Safety

| | |
|---|--|
| Safety coated glassware | Yes, except the evaporating flask |
| Over temperature protection of the bath | Separate monitoring circuit with manual reset Error if temperature is 15 °C above set value |
| Rotation | Soft start Stop in case of blocked rotation |

At any Error

Bath lowering, heater off, rotation off
Type of error showed on display
Reset with main switch

Display

| | |
|------------------------------|--------------|
| Bath temperature | 1 °C steps |
| Cooling temperature (option) | 1 °C steps |
| Vapor temperature | 1 °C steps |
| Set rotation speed | 1 rpm steps |
| Set bath temperature | 1 °C steps |
| Actual vacuum | 1 mbar steps |
| Set vacuum | 1 mbar steps |

Sensors

| | |
|-------------------|---------------------|
| Vapor temperature | PT-1000, 2 wire |
| Bath temperature | PT-1000, 2 wire |
| Vacuum | Ceramic, capacitive |

Features R-250 Pro






| | |
|---|--|
| Two displays | All parameters at a glance on two large displays Always all information available and easy operation at the same time |
| Easy flask handling | Snap flask coupling to fix the evaporating flask Ensures the safe and easy mounting of the evaporating flask by a single person |
| Safety coated glass ware | All glass assemblies (except the evaporating flask) are coated with a robust and transparent safety coating To protect the operator from injuries in case of a glass breakage |
| Indication of process parameters of R-250 Pro | Displays all operating parameters Indicates set and actual values of heating bath, coolant temperature and rotation speed |
| Dynamic distillation | Distillation process starts immediately after choosing the solvent from the library Starts the distillation immediately and adjusts the vacuum dynamically – even if the chiller or bath have not reached their set temperatures. |
| Remote control | The Interface I-300 Pro can be removed and used in combination with a 15 m communication cable to remotely control the complete system Enables the Rotavapor® to be operated from a distance or from behind a closed fume hood |


| | |
|---------------------------|--|
| Remote monitoring | <p>BUCHI Monitor APP for iOS, Android and Windows offers push notifications and live view of all process parameters</p> <p>Allows to track current status of distillation remotely via smart- phones / tablets and informs user when process is terminated</p> |
| Charting | <p>All parameters are graphically displayed to facilitate the supervision of distillation. With the BUCHI Monitor APP also on your mobile device</p> |
| Data recording | <p>Process can be exported on an SD card for further analysis and traceability purpose</p> <p>Enables the continuous recording of all process parameters</p> |
| Wear part library | <p>Internal library lists common wear parts inclusive of order code</p> <p>Allows a convenient replacing process of wear parts and alerts user to check vacuum seal</p> |
| Rotavapor® OpenInterface | <p>Open Interface allows to have an status overview of different BUCHI devices and to have communication between them. Full control on all parameters and functions via PC.</p> |
| Leak test | <p>Integrated test checks system for possible leaks and displays result</p> <p>Allows verifying tightness of the system automatically</p> |
| Different operating modes | <p>Manual vacuum control</p> <ul style="list-style-type: none"> • Manual management of pressure settings and aeration <p>Timer function</p> <ul style="list-style-type: none"> • Manual vacuum control, stops process after preset time has elapsed <p>Continuous pumping</p> <ul style="list-style-type: none"> • Pump runs constantly <p>Methods (SOP's)</p> <ul style="list-style-type: none"> • Performs distillation according to a sequence of programmed steps with defined times and parameters and graphical illustration |
| Multi-languages | en, de, fr, it, es, zh, ja, ru, pt-br, id, ko |
| Overpressure prevention | Automatic aeration when pressure is above 1000 - 1300 mbar (adjustable) |
| ECO-mode | Shuts down activity of bath and chiller, hence lowers energy consumption if system remains inactive for a predefined time period |

Spare parts and accessories





Accessories

System


| | Order no. | Image |
|--|-----------|---|
| <p>Vacuum Pump V-600</p> <p>Chemically resistant 3-stage diaphragm pump. It impresses with its silent and economical operation. Capacity and final vacuum: 3.1 m³/h, 1.5 mbar</p> | 11V600800 |  |
| <p>Vacuum Pump V-600</p> <p>Chemically resistant 3-stage diaphragm pump. It impresses with its silent and economical operation. With secondary condenser. Capacity and final vacuum: 3.1 m³/h, 1.5 mbar</p> | 11V600810 |  |
| <p>Holder vacuum pump</p> <p>To safely store the vacuum pump beneath the bath</p> | 11071091 |  |
| <p>Vacuum valve, 4 mm, 24 V, connection piece 12.5 mm</p> <p>Electrical valve for vacuum regulation when operated with a non-BUCHI vacuum pump.</p> | 11055928 |  |
| <p>Vacuum pump Edwards nXDS10iC cpl.</p> <p>Dry scroll vacuum pump with high chemical stability and a pumping speed of 10 m³/h at 50 Hz and 60 Hz.</p> | 11080656 |  |
| <p>Recirculating Chiller F-360 50 Hz</p> <p>Strong chiller with a cooling capacity of 5600 W at 15°C. Controllable via the Rotavapor R-250 Pro. The temperature range is -10°C to 20 °C and it is operated with 400 V (3PH,N,G) 50 Hz.</p> | 11F36001 |  |



| | Order no. | Image |
|--|-----------|---|
| <p>Recirculating Chiller F-360 60 Hz</p> <p>Strong chiller with a cooling capacity of 6600 W at 15°C. Controllable via the Rotavapor R-250 Pro. The temperature range is -10°C to 20 °C and it is operated with 400 V (3PH,N,G) 60 Hz.</p> | 11F36002 |  |

Safety




| | Order no. | Image |
|--|-----------|---|
| <p>Manual flask handler for 50 L flask</p> <p>For easy mounting and removal of the flasks along with safe transport</p> | 041414 |  |
| <p>Manual flask handler for 20 L flask</p> <p>For easy mounting and removal of the flasks along with safe transport.</p> | 041410 |  |
| <p>Splash protection (cpl.)</p> | 041420 |  |
| <p>Flask crane</p> <p>For the safely secured transport of a 50 liter flask. Incl. the 50 liter manual flask handler.</p> | 041494 |  |

Automation







| | Order no. | Image |
|--|-----------|---|
| <p>Foam detector assembly</p> <p>Internal sensor detects rising foam and triggers a short aeration pulse, eliminating foam.</p> <p>Only in combination with a descending glass assembly.</p> | 11056083 | |
| <p>Cooling water flow sensor</p> <p>Checks the flow of coolant, stopping operation when flow of coolant is insufficient or interrupted.</p> | 11055971 |  |

| | Order no. | Image |
|--|-----------|---|
| <p>Cooling water temperature sensor</p> <p>Needed to display the coolant temperature for optimal distillation settings.</p> | 11055988 |  |
| <p>Level sensor for receiving flask</p> <p>For defined concentration of product or to prevent an overflow of the secondary condenser if combined with a Vacuum Pump V-600 with secondary condenser</p> | 11056192 |  |
| <p>Automatic bath replenishment</p> <p>Detects the water level in the bath to ensures a constant level and protects the bath from running dry.</p> | 11075555 | |

Various


| | Order no. | Image |
|---|-----------|---|
| <p>Vapor duct with integrated sinter plate</p> <p>The integrated sinter plate P3 protects the condenser assembly against powder and dust during the drying process.</p> | 041100 |  |
| <p>Flange adapter for flasks , SJ29.2/32</p> <p>To use a 1, 2 or 3 L evaporating flask with SJ29.2/32</p> | 11058738 |  |
| <p>Stopper, PE, 120 mm</p> <p>To close the evaporating flask</p> | 11057349 | |
| <p>Floating balls</p> <p>Reduces the loss of water in the bath and safes energy.</p> | 035001 | |
| <p>Flask outlet suction system</p> <p>With magnetic tip and PTFE tube (diameter 10 mm)</p> | 041464 | |
| <p>IQ/OQ R-250 Pro</p> <p>official BUCHI document</p> | 11071749 | |
| <p>Communication cable. BUCHI COM, 15 m, 6p</p> <p>Enables connection between Rotavapor®, Interface, Vacuum Pump, Recirculating Chiller, VacuBox and LegacyBox.</p> | 11064090 |  |

Spare parts

| | Order no. | Image |
|---|-----------|---|
| Evaporating flask 50 L D150 For a sample volume of max. 25 L. Flange size 150 mm. | 041339 |  |
| Evaporating flask 20 L D150 For a sample volume of max. 9 L. Flange size 150 mm. | 041432 |  |
| Drying flask 20 L D150 Special shaped to dry up to 9 L product. Flange size 150 mm. | 041393 | |
| Drying flask 50 L D150 Special shaped to dry up to 25 L product. Flange size 150 mm. | 041394 | |
| Vapor duct | 041084 |  |
| Vacuum seal, PTFE base | 041095 |  |
| Evaporating flask seal, complete | 041121 |  |
| FEP coated silicon gasket for DN25 (set of 5) | 11056381 |  |
| FEP coated silicone gasket for DN40 (set of 5) | 11056382 |  |

Performance

The following chart shows the maximum distillation rate.



The maximum achievable distillation rate is not just related to the heating capacity, but also on rotation speed, flask size and temperature difference between bath and cooling.

To get a high distillation rate:

- set a high rotation speed
- choose a large evaporating flask (50 L)
- set a high temperature difference between bath and cooling



11594233C en Rotavapor® R-250 Pro 07.2024
Technical data are subject to change without notice
Quality Systems ISO 9001

We are represented by more than 100 distribution partners worldwide.
Find your local representative at:

www.buchi.com

Quality in your hands
