

## Fat Determination in Meat Products according to Weibull-Stoldt Method *FatExtractor E-500: Soxhlet Extraction after hydrolysis with the HydrolEx H-506*

A simple and reliable procedure for fat determination of meat products according to Weibull-Stoldt is introduced. The sample is hydrolyzed with the HydrolEx H-506. The Soxhlet extraction is performed with the FatExtractor E-500. Calculation of total fat content follows gravimetrically after the extract has been dried to a constant weight. This application follows official methods (eg. EN 98/64/EG, AOAC 963.15, ISO 1443:1973). The combination of the new HydrolEx H-506 and the FatExtractor E-500 increases the sample throughput.

### 1. Introduction

Fat determination is one of the key analysis performed in the food industry. The samples require a hydrolysis step with hydrochloric acid to break the chemically bound and naturally encased fat from the matrix. Afterwards, the fat is extracted with a suitable solvent according to Soxhlet. After the extract has been dried to a constant weight the total fat content is determined gravimetrically.

### 2. Experimental

Equipment: HydrolEx H-506, FatExtractor E-500 Soxhlet

Samples: Cooked sausage LVU No. 16-01j with a certified fat content of 27.46 g/100 g (+/- 0.595 g/100 g); Cervelat sausage, declared fat content 20 g/100g, purchased at a local supermarket; Minced beef meat, declared fat content 13 g/100g, purchased at a local supermarket.

Determination: 20 g of quartz sand was added to a glass sample tube and 2 g Celite® 545 was placed on top. The samples were weighed into a hydrolysis vessel containing 2 g of Celite®. After adding 2 x 50 mL hydrochloric acid (4 M) into each vessel the samples were hydrolyzed for 30 min using the H-506. The hydrolyzate was transferred and the vessels washed with warm (50 °C) deionised water, until a neutral pH was obtained. The glass sample tubes were dried in a vacuum oven, drying oven or microwave oven. After cooling down in a desiccator another layer of quartz sand (20 g) was added to the sample tube. The extraction was performed using the E-500 (Figure 1) applying the parameters specified in Table 1.

Table 1: Parameters for the extraction with the FatExtractor E-500

#### Method parameters

Solvent	Petroleum ether / Diethyl ether
Extraction step	20 cycles (heating level 5 - 6 <sup>1</sup> )
Rinse step	5 min (heating level 5 - 6 <sup>1</sup> )
Smart-Drying	On
Solvent volume	100 mL

The samples were extracted in triplicate. The extracts were dried to a constant weight in a drying oven at 102 °C and the total fat content was calculated.



Figure 1: FatExtractor E-500 SOX

### 3. Results

The determined fat contents are presented in Table 2. The results correspond to the certified values of the reference materials. The determinations show low relative standard deviations.

Table 2: Determined fat content in meat products, fat in g/100 g (relative standard deviation in brackets), n=3

Solvent	Petroleum ether	Diethyl ether
Cooked sausage	27.86 (0.25)	28.06 (0.46)
Cervelat sausage	22.64 (0.45)	22.77 (0.19)
Minced beef meat	12.50 (0.88)	12.60 (0.46)

### 4. Conclusion

The determination of fat in meat products using the HydrolEx H-506 and the FatExtractor E-500 provides reliable and reproducible results. These results correspond well to the labelled values, with low relative standard deviations (rsd).

With the FatExtractor E-500 Soxhlet, the time per cycles is reduced significantly. The programmed 20 cycles are completed in approx. 70 min.

### 5. References

- [1] EN 98/64/EG Commission Directive 98/64/EC Fat in feeding stuffs
- [2] AOAC 963.15 Fat in Cacao Products
- [3] ISO 1443:1973 Meat and meat products -- Determination of total fat content

For more detailed information and safety considerations please refer to the Application Note No. 380/2019.

<sup>1</sup> Heating level proposed by the system depending on the selected solvent.