

Separation of Glycerides by Flash Chromatography Using Different Reversed-Phase Silica Types


Triglycerides are fatty acids present as esters in combination with glycerol. Tristearin, a glyceryl ester of stearic acid, is a simple triglyceride having three identical acyl chains. It is found in both plants and animals and may be used as a drug delivery vehicle for target drug compounds.

This application demonstrates purification of a mixture of glycerides (mono-, di-, and tristearin) using C18 and amino phase chemistries. The C18 uses hydrophobic interaction to determine separation, whereas the amino phase uses lipophilic interaction between the stationary phase and the fatty acyl chain of the analyte to determine separation. The amino phase cartridge can be used in a normal phase or a reversed phase mode, depending on the solvents employed. When used in normal phase mode, the less polar compounds elute first, followed by the elution of the polar one's selectivity.

Experimental



Chromatography system	Pure C-815	C18 cartridge	
Flash cartridge	FlashPure ID C18 12g	Solvent A	Acetonitril
	FlashPure ID Amino 12g	Solvent B	Methylene chloride
Flow rate	30ml/min	Gradient	10-100%: 0-8min
Equilibration	5min	(Solvent B %)	100%: 8-11min
ESLD	Yes	Amino cartridge	
UV wavelengths	254nm, 210nm	Solvent A	Acetonitril
Sample injection	solid	Solvent B	Methylene chloride
		Gradient	100%: 0-2min
		(Solvent B %)	100-10%: 2-11min
			10%: 11-13min

Results



Flash chromatogram of glyceride mixture on a C18 cartridge.

Compound ID:
1. Monostearin
2. Distearin
3. Tristearin



Flash chromatogram of glyceride mixture on an amino cartridge.