

SHORT NOTE

Fat determination in Yoghurt and Curd Cheese using the SpeedExtractor E-916

The work presented is the result of collaboration between BÜCHI and Nestlé. The purpose of this study was to compare PLE results with those gathered by the application of the classical Mojonnier and Weibull-Stoldt (Soxhlet) methods and to skip acid hydrolysis prior to extraction for selected products. The described method has been validated in-house amongst other samples for yoghurt and curd cheese and can be extended to other type of foods. The results of the extraction with the SpeedExtractor E-916 are equivalent to the results achieved by the classical extraction methods (i.e. Mojonnier and Weibull-Stoldt).

Introduction

Pressurized Liquid Extraction (PLE) is a well-established alternative for extracting lipids and/or various components faster and cheaper, but less frequently used for the verification of fat labels in food products. A reliable and fast procedure using PLE to determine the fat content in samples with easily accessible fat and high moisture content (above 50 %) is presented below. The extraction was carried out using the SpeedExtractor E-916 in combination with MultivaporTM P-12 for solvent evaporation in parallel. The fat content was determined gravimetrically.

Experimental

Instrumentation: SpeedExtractor E-916, Multivapor™ P-12 with Vacuum pump V-700 and Controller V-855, drying oven

HIRZ Rahnquark

Figure 1: Curd cheese

Samples: Yoghurt Nestlé LC1 nature, curd cheese

The samples were ground with poly(acrylic acid) and diatomaceous earth in a mortar to a homogeneous powder and quantitatively transferred into a 20 ml extraction cell. The extraction was carried out with the SpeedExtractor E-916 using the parameters shown in Table 1. The samples were extracted in triplicate. Per position (sample) about 30 ml of solvent was used.

Table 1: Extraction method for fat determination with SpeedExtractor E-916

Table 1. Extraction method for lat determination with SpeedExtractor E-910		
Temperature	100 °C	
Pressure	100 bar	
Solvent	n-Hexane 90 %; Ethanol 10 %	
Cells	20 ml	
Vials	60 ml	
Cycles	2	
Heat-up	1/1 min	
Hold	10/10 min	
Discharge	2/2 min	
Flush with solvent	1 min	
Flush with gas	3 min	
Total time	42 min	

During the extraction the extracts of each sample were collected simultaneously in 60 ml vials. The solvent was evaporated to dryness using the MultivaporTM P-12 with the pressure gradient shown in Fig 2.




Figure 2: Pressure gradient for the evaporation of the solvent using Multivapor™ P-12

The extracts were dried to a constant weight in a drying oven (102 °C) and the fat contents were calculated.

Results

The fat content of yoghurt and curd cheese are shown in Table 2. The results correspond to the reference values.

Table 2: Determined fat contents in food samples, fat in g/100 g (relative standard deviation for the experimental results in brackets). n = 3

Sample	Expected fat contents	Experimental results
Yoghurt	3.25	3.18 (1.03 %)
Curd cheese	15.07	14.96 (0.50 %)

Conclusion

The determination of the fat content in these dairy products by PLE using SpeedExtractor E-916 provides reliable and reproducible results which correspond to the expected values. This method increases the sample throughput and is more environmental friendly due to less solvent consumption compared to Soxhlet extraction.

Acknowledgement

We sincerely thank Nestlé Research Center, Switzerland for their support in developing this application note.

References

SpeedExtractor E-916 operation manual

For more detailed information refer to Application note 089/2012