




## SyncorePlus

### Pre-installation checklist

The following checklist will assist you in site preparation for your BUCHI instrument. Careful attention to the requirements presented in this document will help simplify instrument installation and ensure that the instrument operates as it was designed.





# **SyncorePlus**

## Pre-installation checklist

All appropriate items in this checklist must be completed, signed and returned to BUCHI. Return your checklist as soon as your laboratory is ready and functional for prompt installation of your SyncorePlus. If you are performing self-installation, please check the box for self-installation on page 4 and return the form to BUCHI.

### Environmental requirements

Place the SyncorePlus in a well ventilated room or fume hood.

Ambient temperature of 5°C to 40 °C

Relative humidity  $\leq$  70%

### Space requirements

Enough space as defined below: (Width x Depth x Height)

For the SyncorePlus: 500 x 520 x 325 mm.

For the optional Recirculating Chiller F-314: 400 x 500 x 660 mm.

For other BUCHI Recirculating Chillers; consult the relevant technical data sheet.

For the optional Vacuum Pump V-300 & Interface I-300: 200 x 325 x 315 mm.

NB Include additional depth for optional secondary condensers.

Keep the recirculating chiller as close as possible to the SyncorePlus to ensure that the coolant tubing can be kept to a minimum length.

Allow at least 40 cm of space between the recirculating chiller and the wall and 40 cm in front of the chiller for sufficient airflow and proper cooling of the chiller.

Allow additional space around the instrument for a comfortable working environment, preparing samples, removing glassware, etc. Take into account any additional obstacles in the area such as gas taps, water taps, drains, sloping ceilings in enclosures, etc...



# **SyncorePlus**

## Pre-installation checklist

### Location of the unit

The unit must be installed on a clean, stable, and level surface.

If using a mains water supply for cooling, then the system should also be within 1.0m of the cooling water supply. Cooling water needs to operate with unobstructed flow with no more than 2 bar of pressure. The cooling water tap needs to have a connection suitable for pushing on and securing a 6mm ID / 9mm OD hose. A drainage sink able to keep secure a drainage hose should also be available.

### Power requirements

A SyncorePlus requires a standard 230V mains electricity socket.

Each accessory (Vacuum Pump, Recirculating Chiller) requires its own standard 230V mains electricity socket.

All mains electricity cables and plugs are supplied.

### Other requirements

When using a Recirculating Chiller, take into account the chiller temperature setting when selecting a coolant. A 40:60 ethylene glycol to water mixture that does not freeze until -16°C is recommended. Allow 6L for an F-305, 8L for an F-308 or 10L for an F-314; this allows enough to fill the chiller and condenser while also leaving enough for future use.

### Facility sign-off

Return your checklist as soon as your facility is ready and functional for prompt installation of your SyncorePlus

Environmental Requirements as specified for your facility

Space Requirements

Other Requirements

Power Requirements

Self-Installation



## **SyncorePlus**

### Pre-installation checklist

Please enter the following information and sign below:

Company name

Customer name

Contact number

Purchase order no.

I certify that the preparations required for the installation of the SyncorePlus been completed and functional as described in this document.

**I further acknowledge that if any of these requirements are not ready by the agreed upon installation date and I have not rescheduled the installation, then I will be responsible for any additional costs relating to the installation of the system.**

Customer signature

Date

Please sign and return by fax or email to BUCHI as soon as your laboratory is ready for installation:

**Attention:**

Find your BUCHI affiliate or partner on:

[www.buchi.com/en/worldwide](http://www.buchi.com/en/worldwide)

---

If you have any questions or queries, then please do not hesitate to contact us.

---

## BUCHI Affiliates:

### Europe

<b>Switzerland/Austria</b> <b>BÜCHI Labortechnik AG</b> CH – 9230 Flawil T +41 71 394 63 63 F +41 71 394 64 64 buchi@buchi.com www.buchi.com	<b>Benelux</b> <b>BÜCHI Labortechnik GmbH</b> Branch Office Benelux NL – 3342 GT Hendrik-Ido-Ambacht T +31 78 684 94 29 F +31 78 684 94 30 benelux@buchi.com www.buchi.com/bx-en	<b>France</b> <b>BUCHI Sarl</b> FR – 91140 Villebon-sur-Yvette T +33 1 56 70 62 50 F +33 1 46 86 00 31 france@buchi.com www.buchi.com/fr-fr	<b>Germany</b> <b>BÜCHI Labortechnik GmbH</b> DE – 45127 Essen T +800 414 0 414 0 (Toll Free) T +49 201 747 49 0 F +49 201 747 49 20 deutschland@buchi.com www.buchi.com/de-de
<b>Italy</b> <b>BUCHI Italia s.r.l.</b> IT – 20010 Cornaredo (MI) T +39 02 824 50 11 F +39 02 575 12 855 italia@buchi.com www.buchi.com/it-it	<b>Russia</b> <b>BUCHI Russia/CIS</b> RU – 127006 Moscow T +7 495 36 36 495 russia@buchi.com www.buchi.com/ru-ru	<b>Spain</b> <b>BUCHI Ibérica S.L.U.</b> ES – 08960 Barcelona T +34 936 06 8010 iberica@buchi.com www.buchi.com/es-es	<b>United Kingdom</b> <b>BUCHI UK Ltd.</b> GB – Suffolk CB8 7SQ T +44 161 633 1000 F +44 161 633 1007 uk@buchi.com www.buchi.com/gb-en

### NIR-Online

**Germany**  
**NIR-Online GmbH**  
 DE – 69190 Walldorf  
 T +49 6227 73 26 60  
 F +49 6227 73 26 70  
 nir-online@buchi.com  
 www.nir-online.de

### America

<b>Brazil</b> <b>BUCHI Brasil Ltda.</b> BR – Valinhos SP 13271-200 T +55 19 3849 1201 F +55 19 3849 2907 brasil@buchi.com www.buchi.com/br-pt	<b>USA/Canada</b> <b>BUCHI Corporation</b> US – New Castle, DE 19720 T +1 877 692 8244 (Toll Free) T +1 302 652 3000 F +1 302 652 8777 us-sales@buchi.com www.buchi.com/us-en
---	--

### Asia

<b>China</b> <b>BUCHI China</b> CN – 200233 Shanghai T +86 21 6280 3366 F +86 21 5230 8821 china@buchi.com www.buchi.com/cn-zh	<b>India</b> <b>BUCHI India Private Ltd.</b> IN – Mumbai 400 055 T +91 22 667 75400 F +91 22 667 18986 india@buchi.com www.buchi.com/in-en	<b>Indonesia</b> <b>PT. BUCHI Indonesia</b> ID – Tangerang 15332 T +62 21 2222 9478 F +62 21 5568 7906 indonesia@buchi.com www.buchi.com/id-in	<b>Japan</b> <b>Nihon BUCHI K.K.</b> JP – Tokyo 110-0008 T +81 3 3821 4777 F +81 3 3821 4555 nihon@buchi.com www.buchi.com/jp-ja
<b>Korea</b> <b>BUCHI Korea Inc.</b> KR – Seoul 153-782 T +82 2 6718 7500 F +82 2 6718 7599 korea@buchi.com www.buchi.com/kr-ko	<b>Malaysia</b> <b>BUCHI Malaysia Sdn. Bhd.</b> MY – 47301 Petaling Jaya, Selangor T +60 3 7832 0310 F +60 3 7832 0309 malaysia@buchi.com www.buchi.com/my-en	<b>Singapore</b> <b>BUCHI Singapore Pte. Ltd.</b> SG – Singapore 609919 T +65 6565 1175 F +65 6566 7047 singapore@buchi.com www.buchi.com/sg-en	<b>Thailand</b> <b>BUCHI (Thailand) Ltd.</b> TH – Bangkok 10600 T +66 2 862 08 51 F +66 2 862 08 54 thailand@buchi.com www.buchi.com/th-th

### BUCHI Support Centers:

<b>South East Asia</b> <b>BUCHI (Thailand) Ltd.</b> TH-Bangkok 10600 T +66 2 862 08 51 F +66 2 862 08 54 bacc@buchi.com www.buchi.com/th-th	<b>Middle East</b> <b>BÜCHI Labortechnik AG</b> UAE – Dubai T +971 4 313 2860 F +971 4 313 2861 middleeast@buchi.com www.buchi.com	<b>Latin America</b> <b>BUCHI Latinoamérica S. de R.L. de C.V.</b> MX – Mexico City T +52 55 9001 5386 latinoamerica@buchi.com www.buchi.com/es-es
---	--	---

We are represented by more than 100 distribution partners worldwide.  
 Find your local representative at: [www.buchi.com](http://www.buchi.com)