



BUCHI Solutions

# Natural Products



## Workflow Natural Products

Natural products have served as a major source for pharmaceutical drugs since many years. Besides that, they play a key role as ingredients for food and cosmetic applications. Due to their enormous variety, the methods used for development and production of these compounds is manifold.



**Discovery**

The discovery of a natural product includes several steps. The process begins with production of a crude extract from a natural source, which is used for a biological screening to identify a bioactive compound. Then, any target of interest gets further processed to achieve a single compound of high purity, which is tested thoroughly for desirable characteristics and functions.



**Development**

Once an active ingredient has been identified, the development of the production process can begin. Here, the process is optimized on a larger scale to ensure a trouble-free production process. Before being incorporated into a final formulation, ingredients need to be preformulated in either a liquid or solid form to facilitate storage. Formulation is used to protect the target material as well as for controlled or target release.



**Testing**

Natural products must undergo strict application testing before they enter the production process. The final product is typically tested for chemical stability and physical integrity under specific storage, transport, and use conditions. In case of pharmaceutical application, clinical trials can take several years to complete and involve evaluating the drug mixture with animal, in vitro, and in vivo studies. If the product is proven safe and efficient, it can be released for production.



**Production**

Inspection of incoming goods and final quality control ensures the desired quality of any natural product. Additionally, quality control during the production helps to analyze the identity, purity, and content of bulk material, intermediates, impurities, and degradation products. Close monitoring of the process enables safe operation and ensures the product meets the required specifications.

# Natural Products R&D Discovery

Extraction

Concentration

Separation

Drying

Analysis

**Solvent Extraction / Soxhlet**

**Evaporation**

**Flash Chromatography,  
Prep HPLC, Prep SFC**

**Freeze Drying**

**Melting Point**



UniversalExtractor  
E-800

SpeedExtractor E-914 /  
E-916

Rotavapor®

Rotavapor®

SyncorePlus

Pure & Pure Essential,  
Sepiatec SFC, Consumables

Lyovapor™

Melting Point

Extraction is used to separate the sample matrix (plant raw material) from the desired compounds. Based on the properties of the target analyte, the solvent choice, the applied pressure, temperature and extraction parameters are essential to reach high extraction efficiency.

Versatile and highly automated extraction instruments are used for the screening process for the evaluation of the raw materials. Soxhlet extraction can be performed also on a Rotavapor® when using the adapted glassware, simplifying the glassware set up and process.

Since extraction requires large amounts of solvent, a concentration step is required prior to downstream processing. Here, rotary evaporation is used to remove the solvent and concentrate the compound of interest.

The use of parallel evaporation can speed up the concentration of multiple samples. Many samples are evaporated simultaneously, which increases sample throughput.

Flash and prep HPLC (preparative high pressure liquid chromatography) & prep SFC (preparative super critical fluid chromatography) are commonly used to purify target compounds: flash is used as a pre-purification step, whereas prep HPLC and SFC increase the purity of the target compound to the maximum.

Following the separation process, molecules of interest are highly diluted and must be concentrated or dried before proceeding with the next steps. Freeze drying can be used to remove solvent from heat-sensitive products with minimal damage.

Melting point analysis can be used to perform quality control on the compound of interest. Determination of the novel compound's melting point serves as a useful indicator of the material's purity.

Application

Features

- Five different extraction methods, including true Soxhlet
- Simultaneous operation of different extraction methods
- Visible extraction process
- Gentle drying of extract

- Extraction under elevated temperature and pressure
- Very efficient extraction as the pressure breaks up the cell structure
- Gentle extraction possible due to adjustment of temperature below the boiling point

- Soxhlet accessory for Soxhlet extraction
- One instrument fits several application: extraction and evaporation
- Extraction chamber size: 200 mL and 500 mL

- Evaporation of a single sample with volumes of 50 - 5000 mL
- Fully communicating system to avoid downtime: solvent library, dynamic distillation, drying mode, leak test, foam sensor
- Foam sensor accessory to automatically handle foam
- Dewar accessory for freeze drying sample preparation


- Multiple samples in the range of 0.5 - 500 mL can be concentrated or dried simultaneously
- Flushback Module to achieve highest analyte recovery and most reliable results
- Interchangeable racks and volume versatility

- Flash instruments for basic or advanced applications
- Flash and prep HPLC in one system (optionally)
- Integrated UV and ELS detection (optionally)
- Compact prep SFC instruments for any sample size
- Compatible with a wide range of flash cartridges, prep HPLC & SFC columns and glass columns

- High-quality freeze drying of aqueous samples (-55 °C, 6 kg)
- Easy way of controlling and monitoring of the freeze drying process

- Automatic determination of melting and boiling points
- Compliant with Pharmacopeia methods (European, USP, Japanese)
- Observation and replay of the phase transition using color display and video recording
- Parallel measurement of up to 3 samples

# Natural Products R&D Development



# Natural Products Production

Sample Preparation for Analysis


Analysis

Analysis

## Extraction

## Melting Point

## NIR Spectroscopy



UniversalExtractor E-800

SpeedExtractor E-914 / E-916

Melting Point

ProxiMate™

NIRFlex N-500

NIR-Online Process Analyzer

Application

Extraction is used as sample preparation for the quality control of final phyto-pharmaceuticals. The quantification of the compound of interest and the determination of possible contaminants or undesired residues are vital to ensure the safety of the final product as well as of the raw material.

Melting point analysis can be used to perform quality control on the compound of interest. Determination of the novel compound's melting point serves as a useful indicator of the material's purity.

NIR is a fast and safe screening method that provides reliable results within (milli)seconds. The non-destructive technique is applied for raw material identification during incoming inspection as well as for quantification of multiple parameters in intermediate and final products. With NIR Atline and Online solutions, you can ensure your quality, optimize all process steps and thus improve your productivity and profitability. The laboratory device can moreover be applied to check content uniformity in tablets and capsules.

Features

- Fully automated solvent extraction
- Suitable for a wide range of organic solvent or water
- High sample throughput by extracting up to 6 samples in parallel

- Speed and throughput combined: able to extract 6 samples in 20 mins
- Reduced solvent consumption gives low running costs
- Complementary workflow with parallel evaporation. Compatible glassware makes extract transfer obsolete.

- Automatic determination of melting and boiling points
- Compliant with Pharmacopeia methods (European, USP, Japanese)
- Observation and replay of the phase transition using a color display and video recording
- Parallel measurement of up to 3 samples

- Integrated touchscreen, stand-alone
- IP69 certified, At-Line instrument
- AutoCal for calibration development without expert knowledge
- Up- and Down-View possible depending on the individual requirements
- Colour measurements according to CIE norms

- Modular design that meets individual requirements
- Compliant with 21 CFR Part 11 as well as US, EU, and Japanese pharmacopoeias regulations
- Highest spectrometer precision
- Identification of raw materials
- Quantification of raw materials, intermediates, and final products

- Diode array technology for high-speed measurement of fast-moving goods
- Industry proven, robust design without moving parts copes with rough conditions e.g., vibrations, extreme temperatures, wind or humidity
- Certified for utilization in potentially explosive gas and dust atmospheres



# Pressurized Solvent Extraction


## Product Details [↗](#)



### System Portfolio & Technical Features


SpeedExtractor	Number of parallel samples <a href="#">↗</a>	Sample volume [mL]	Number of solvents	Collection vials
E-916	6	10, 20, 40	2 or 4	60 mL, 220 mL, 240 mL, SyncorePlus R-12, Waste
E-916 XL	6	60	2	60 mL, 220 mL, 240 mL, SyncorePlus R-12, Waste
E-914	4	10, 20, 40, 80, 120	2 or 4	60 mL, 220 mL, 240 mL, SyncorePlus R-12, SyncorePlus R-6, flasks for Rotavapor®, Waste

### Accessories




#### Extraction Cells

Adjust the cell size to your sample amount. Extraction cells are available for 10, 20, 40, 60, 80 and 120 mL.




#### Collection Vessels

Full compatibility of SpeedExtractor collection glassware with SyncorePlus Analyst and Multivapor™.




#### Carrier for Extraction Cells

Used for sample preparation and carrying of all samples ensuring a safe handling.



#### SpeedExtractor Record Software

Allows for managing and storage of methods, graphical visualization and documentation.



#### IQOQ Documentation

Installation / Operational Qualifications documentation.

### Automation and Workflow



#### SyncorePlus Analyst

Parallel evaporation for pre-analytical and fast concentration achieving high sample throughput.



#### Multivapor™ P-6 / P-12

Efficient evaporation for multiple samples. With its ease of use, the Multivapor™ is tailored to maximize efficiency.



# Solvent Extraction


## Product Details [↗](#)



### System Portfolio & Technical Features


UniversalExtractor Model	Extraction method <a href="#">▶</a>	Chamber heater	Inert gas supply	Large sample volume
E-800 Standard	Soxhlet, Continuous flow	no	no	yes (option)
E-800 Pro	Soxhlet, Soxhlet warm, Hot extraction Continuous flow, Twisselmann	yes	no	yes (option)
E-800 Pro Inert	Soxhlet, Soxhlet warm, Hot extraction Continuous flow, Twisselmann	yes	yes	yes (option)

### Accessories and Consumables




#### Holders and Support

Beneficial holders facilitate the easy handling of the beakers and vessels.




#### LSV Glass Parts

Expand your sample volume by using the Large Sample Volume option (LSV).




#### Extraction Thimbles

Best quality and optimized dimensions for sample extraction.



#### Extraction Reports App

Monitoring and reporting of your extraction process.



#### IQOQ Documentation

Installation / Operational Qualifications documentation.

### Automation and Workflow



#### Recirculating Chiller F-305, F-308, F-314

Convenient central temperature setting, energy-saving ECOmode, and automatic start / stop.



#### Mixer B-400

Optimized sample preparation, combines homogenization efficiency with simple operation.



#### SyncorePlus Analyst

Parallel evaporation for pre-analytical and fast concentration for high sample throughput.



# Laboratory Evaporation


## Product Details [↗](#)



### System Portfolio & Technical Features

Rotavapor® Model	Max. Flask Size	Lift System	Vacuum Pump	Chiller
R-300 <a href="#">↗</a>	1 or 5 L	electric or manual	no	no
Rotavapor® system RS-300	1 or 5 L	electric or manual	yes (final vacuum: 5 mbar)	yes (optional)
R-100 <a href="#">↗</a>	4 L	manual	no	no
Rotavapor® system RS-100	4 L	manual	yes (final vacuum: 10 mbar)	yes (optional)

### Accessories




**Extraction glassware**  
To perform Soxhlet extraction on the Rotavapor. Available extraction chambers: 200 and 500 mL.



**Beaker Flasks**  
Beaker flasks with large screw-cap opening for easy retrieval of substances.



**Foam Sensor**  
Prevents sample from foaming into the condenser by automatic short aeration of the system.



**Monitor App**  
Allows to monitor all BUCHI instruments on a mobile device thanks to the BUCHI Bluetooth Dongle.

### Related Instruments



**Vacuum pumps V-300, V-600**  
Chemically resistant diaphragm pump. Final vacuum: V-300: 5 mbar / V-600: 1.5 mbar.



**Recirculating Chiller F-305, F-308, F-314**  
Convenient central temperature setting, energy-saving ECOMode, and automatic start / stop.

### Condensers



A	C	V	S	CR	E	BY	HP
Diagonal	Cold trap	Vertical	Reflux	Cold trap reflux	Expansion	Double jacket	High performance
•	• •	• •	•	• •	•	•	• •

- Reduced height
- Low boiling point
- Standard
- Foaming products
- Reflux reactions
- Increased distillation rate

# Parallel Evaporation


## Product Details [↗](#)



### System Portfolio & Technical Features


SyncorePlus Model	Method	Number of Samples	Max Temperature	Rotation
Analyst <a href="#">↗</a>	concentration to final volume	4, 6, 12	100 °C	60 – 400 rpm
Polyvap <a href="#">↗</a>	evaporation to dryness	4, 6, 12, 24, 48, 96	100 °C	60 – 400 rpm

### Accessories




#### Graduated Vessels

Concentrates of up to 12 samples to a predefined residual volume using SyncorePlus Analyst, ranging from 0.3 mL, 1 mL, to 3 mL per sample.




#### Interchangeable Racks

Broad selection of racks for SyncorePlus Polyvap to accommodate any workflow or throughput; that come in 4, 6, 12, 24, 48 and 96 positions.




#### Flushback Module for R-6, R-12

Improves recovery rates by flushing back analyte that adheres to the glass walls with condensed vapor during the evaporation process.



#### Solid Phase Extraction Cover

All essential work-up steps including evaporation of the eluates are achieved without any sample handling between the individual steps.



#### IQOQ Documentation

Installation / Operational Qualifications documentation.

### Automation & Workflow



#### Interface I-300 Pro

Full automation with the Interface I-300 Pro. Start the process and walk away. Specified evaporation methods and predefined specific solvents.



#### Recirculating Chiller F-305, F-308, F-314

Convenient central temperature setting, energy-saving ECOMode, and automatic start / stop.



#### Vacuum Pump V-300

Essential system component for quiet, eco-friendly operations and the ability to distill solvents of any volume and with any boiling point.



#### SpeedExtractor E-914 / E-916


Pressurized Solvent Extraction (PSE) for increased productivity by processing up to 6 samples in parallel. Streamlined workflow of the sample preparation thanks to ease of sample loading and ready to use extract collection.

# Industrial Evaporation


## Product Details [↗](#)




### Accessories




**Extraction glassware**  
To perform cold extraction on the Rotavapor R-220 Pro. Available extraction chamber: 4 L.



**Manual Flask Handler**  
For easy mounting and removal of the 20 L and 50 L flasks along with safe transport.



**Foam Sensor**  
Detects rising foam and triggers a short aeration pulse, eliminating foam.



**Monitor App**  
Allows to monitor all BUCHI instruments on a mobile device thanks to the BUCHI Bluetooth Dongle.

### Related Instruments



**Vacuum Pump V-600**  
Chemically resistant 3-stage diaphragm pump.  
Final vacuum: 1.5 mbar




**Recirculating Chiller F-325**  
Chiller that also operates as a trolley and host of the Vacuum Pump V-600 for the R-220 Pro.

### Condensers

R-220 Pro

R-250 Pro



D D2 DB2 DB RB R C D2 D3 RB2 R2



Height 175 cm 175 cm 150 cm 150 cm 143 cm 163 cm 158 cm 230 cm 230 cm 210 cm 226 cm

- Low boiling points and / or foaming products
- Minimum emissions
- High boiling points
- Reflex reactions
- Very low boiling point
- Reduced height

### System Portfolio & Technical Features

Rotavapor® Model	Sample Size (per batch)	Heating Bath Temperature	Distillation Rate of Ethanol	EX Protection
R-220 Pro <a href="#">↗</a>	max 12 L	up to 180 °C	up to 18 L / h	no
R-250 Pro <a href="#">↗</a>	max 30 L	up to 180 °C	up to 19 L / h	no
R-220 EX <a href="#">↗</a>	max 12 L	up to 150 °C (T3)	up to 12 L / h	yes
R-250 EX <a href="#">↗</a>	max 30 L	up to 150 °C (T3)	up to 19 L / h	yes



# Flash Chromatography / Prep HPLC


## Product Details [↗](#)



### System Portfolio & Technical Features

	Pure C-805 Flash <a href="#">↗</a>	Pure C-810 Flash	Pure C-815 Flash	Pure C-830 Prep	Pure C-835 Prep	Pure C-850 FlashPrep
Mode	Flash	Flash	Flash	prep HPLC	prep HPLC	Flash & prep HPLC
Flow Rate (flash mode)	250 mL / min	250 mL / min	250 mL / min			250 mL / min
Flow Rate (prep HPLC mode)				100 mL / min	100 mL / min	100 mL / min
Max Pressure (flash mode)	50 bar	50 bar	50 bar			50 bar
Max Pressure (prep HPLC mode)				300 bar	300 bar	300 bar
UV Scan Function <a href="#">↗</a>	no	yes	yes	yes	yes	yes
UV-Vis Wavelengths Range	200 – 400 nm	200 – 800 nm	200 – 800 nm	200 – 800 nm	200 – 800 nm	200 – 800 nm
ELSD <a href="#">↗</a>	no	no	yes	no	yes	yes

### Accessories




#### Dry Air Supply

Production of particle-free dry air on demand without the need for operator attention. Air is used in the Pure system to carry sample to the ELSD and to purge cartridges and sample loaders.




#### Solid Loader

Easy and flexible filling of the sample (max. 40 g) in an empty tube and connection directly to the Pure system. This equipment can handle up to 50 bar (725 psi) and therefore provides maximum flexibility.




#### Racks

Several different sizes of racks and glass tubes, which enable an optimal collection depending on the fractions sizes. All racks get automatically identified by the Pure system via an RFID tag.




#### XL Consumable Holders

Allows the connection of big size cartridges (750 – 5000 g) and prep HPLC columns (50 – 70 mm ID) with Pure.




#### Solvent Platform

Extra solvent platform on top of the Pure system which provides space for four bottles. This allows for better use of available space and reduces risk of spills.



#### Sampling Pump

For large sample volumes injected manually on a large flash cartridge. Flow rates up to 250 mL / min and 50 bar.



#### IQOQ Documentation

Installation / Operational Qualifications documentation.

### Consumables



#### FlashPure

A wide range of flash cartridges

Phases	Silica, C18, Amino, Diol, C18 WP, Alumina
Particle sizes	15 – 50 µm
Particle shapes	Irregular, spherical
Cartridge sizes	4 – 5000 g



#### PrepPure

Highest performance for prep HPLC applications

Phases	Silica, C18, C18WP, C4WP, C18AQ, Diol, PEI, 2-EP, Chiral phases
Particle sizes	5 – 15 µm
Particle shapes	Spherical
Column sizes	4.6 – 70 mm ID, 150 & 250 mm Lengths



#### GlasPure

Scale-up purification

Lengths	100 – 900 mm
IDs	15 – 100 mm
Silica capacities	9 – 3400 g



# Flash Chromatography (basic & modular)

Product Details [↗](#)




## System Portfolio & Technical Features

	Pure Chromatography C-900 <a href="#">↗</a>	Pure Chromatography C-900 + Fraction Collector <a href="#">↗</a>	Pure Chromatography C-900 + UV Detector <a href="#">↗</a>	Pure Essential Chromatography System <a href="#">↗</a>
Modules	<ul style="list-style-type: none"> <li>3 piston flash pump with binary gradient and a controller with novel and state-of-the-art software</li> </ul>	<ul style="list-style-type: none"> <li>3 piston flash pump with binary gradient and a controller with novel and state-of-the-art software</li> <li>Enclosed fraction collector with active ventilation</li> </ul>	<ul style="list-style-type: none"> <li>3 piston flash pump with binary gradient and a controller with novel and state-of-the-art software</li> <li>Highly sensitive 4 times fixed wavelengths UV detector</li> </ul>	<ul style="list-style-type: none"> <li>3 piston flash pump with binary gradient and a controller with novel and state-of-the-art software</li> <li>Enclosed fraction collector with active ventilation</li> <li>Highly sensitive 4 times fixed wavelengths UV detector</li> </ul>

### Pure Essential Chromatography System


Mode	Flash
Flow rate	300 ml / min
Max pressure	50 bar
UV wavelengths	4
UV spectra	254, 275, 325, 365 nm

## Accessories




### Racks

The Pure fraction collector can be equipped with two racks and is compatible with several different-sized glass tubes, funnels, and flasks, enabling an optimal collection depending on the fraction size.




### Solid loader hardware

The Pure Solid Loader allows for easy and flexible filling of the sample in an empty tube and connection directly to the Pure Essential Chromatography system. This equipment can handle up to 50 bar (725 psi), providing maximum flexibility.




### Solvent platform

The Pure fraction collector can be equipped with an extra solvent platform on top which provides space for four bottles. This allows for better use of available space and reduces risk of spills.




### Loops and chambers

A liquid sample can be injected manually on the flash cartridge or loaded into a loop or a chamber, depending on the sample size.




### Mixing chambers

The Pure high-pressure mixing chambers ensure proper mixing of the solvents, especially when using low concentrations.



### Cartridge holder

The Pure Essential Chromatography system can handle different flash cartridges (sizes vary between 4 and 5000 g) or glass columns (ID up to 100 mm and lengths up to 900 mm) depending on the quantity of sample purified.



### IQOQ Documentation

Installation / Operational Qualifications documentation.

## Consumables



### FlashPure

A wide range of flash cartridges

Phases	Silica, C18, Amino, Diol
Particle sizes	20 – 50 µm
Particle shapes	Irregular, spherical
Cartridge sizes	4 – 5000 g



### GlasPure

Scale-up purification

Lengths	100 – 900 mm
IDs	15 – 100 mm
Silica capacities	9 – 3400 g



## Prep SFC

### Product Details [↗](#)



### System Portfolio & Technical Features

	Sepiatec SFC-50 <a href="#">↗</a>	Sepiatec SFC-250 <a href="#">↗</a>	Sepiatec SFC-660 <a href="#">↗</a>
Mode	prep SFC	prep SFC	prep SFC
<b>Columns</b>			
Dimensions	4 – 16 mm ID, max 250 mm lengths	15 – 30 mm ID, max 250 mm lengths	30 – 50 mm ID, max 800 mm lengths
<b>Pumps</b>			
CO <sub>2</sub> pump flow rate and pressure	Max 30 mL / min, 400 bar	Max 150 mL / min, 400 bar	Max 400 mL / min, 400 bar
Max operating pressure	400 bar	400 bar	400 bar
Total flow rate At 40 % modifier	50 mL / min	250 mL / min	660 mL / min
<b>Detectors</b>			
DAD wavelengths range	190 – 720 nm	190 – 720 nm	190 – 720 nm
DAD wavelengths selectable	8	8	8
<b>Sample injection</b>			
Syringe	1 mL Other volumes on request	5 mL Other volumes on request	25 mL Other volumes on request
Loop	0.5 mL	2.5 mL	10 mL
Stack injection	Standard	Standard	Standard
<b>System controller</b>			
Back pressure regulator	80 – 250 bar	80 – 250 bar	80 – 250 bar

### Accessories

#### Add-on pump

The add-on pump is useful for separations at modifier concentrations below 10 %. The pump adds modifier to the separated sample and avoids precipitation of the sample in the gas-liquid separator.

#### Chiller

The chiller cools the pump heads of the CO<sub>2</sub> pump and is connected simultaneously to the pre-cooling and the CO<sub>2</sub> pump. The cooling is needed to keep the CO<sub>2</sub> in a liquid state.

#### Mass spectrometer (MS)

All Sepiatec SFC instruments can be connected with an external MS detector. Further info are given on request.

#### Evaporative light scattering detector (ELSD)

All Sepiatec SFC instruments can be connected with an external ELS detector. Further info are given on request.

### Consumables



#### PrepPure

Highest performance  
for prep HPLC & SFC  
applications

#### Phases

Silica, C18, C18WP, C18AQ, C4WP,  
Diol, 2-EP, PEI, Chiral phases

#### Particle sizes

5 – 15 µm

#### Particle shapes

Spherical

#### Column sizes

4.6 – 70 mm ID, 150 & 250 mm Lengths

### Chiral phases

#### Phase

##### Immobilized polysaccharides

iADMPC (Amylose tris-(3,5-dimethylphenylcarbamate))

iCDMPC (Cellulose tris-(3,5-dimethylphenylcarbamate))

iCDCPC (Cellulose tris-(3,5-dichlorophenylcarbamate))

##### Coated polysaccharides

cCDMPC (Cellulose tris-(3,5-dimethylphenylcarbamate))

cADMPC (Amylose tris-(3,5-dimethylphenylcarbamate))

##### Brush type

iBT (immobilized brush-type phase)




# Freeze Drying


## Product Details [↗](#)




### Accessories [↗](#)




**Glassware**  
Large variety of glassware such as beaker flasks to meet laboratory needs.



**Lyovapor Software**  
Easy way of controlling and monitoring the Freeze Drying process and generating reports.



**Edwards Scroll Pump**  
Stable and high-quality vacuum pump for freeze drying of organic solvents.



**Sensors**  
Choice of product temperature and vacuum sensors for end point determination.

### Drying Chambers

The Lyovapor™ has a large choice of drying chambers that can be used in any combination.




Acrylic chamber with heated shelves and stoppering for vial drying.



Acrylic chamber with 12 manifold top and shelves for bulk, vial and flask drying.



Acrylic chamber with heated shelves for bulk and vial drying.




Manifold with 12, 24 or 36 valves for flask or beaker drying.

### System Portfolio & Technical Features

Lyovapor Model	Condenser Capacity	Lowest Condenser Temperature	Condensing Capacity	Min. System Vacuum	Drying Shelf Temperature
L-200 <a href="#">↗</a>	6 kg	-55 °C + / -3 °C	6 kg / 24 h	30 mTorr / 0.04 mbar	Up to 60 °C + / -1 °C
L-250 <a href="#">↗</a>	5 kg	-85 °C + / - 2 °C	4 kg / 24 h	30 mTorr / 0.04 mbar	Up to 60 °C + / -1 °C
L-300 <a href="#">↗</a>	Infinite	-105 °C + / -3 °C	12 kg / 24 h	30 mTorr / 0.04 mbar	Up to 60 °C + / -1 °C


# Spray Drying Product Details



## System Portfolio & Technical Features


	Nano Spray Dryer B-90 HP 	Mini Spray Dryer S-300 
Particle Size	200 nm – 5 µm	2 – 25 µm (60 µm with ultrasonic package)
Particle Nature	dry	dry
Particle Size Distribution	narrow	broad
Max. Sample Throughput	200 mL / h	1 L / h
Min. Sample Amount	200 mg / 2 mL	5 g / 10 mL
Yield	up to 90 %	up to 70 %
Sample Viscosity	up to 5 cps	up to 300 cps
Sample Composition	aqueous and organic solutions, suspensions or emulsions NO acidic or alkaline	aqueous, organic, acidic and alkaline solutions, suspensions or emulsions NO acidic or alkaline organic mixture
Application	drying micronisation, agglomeration matrix encapsulation amorphous solid dispersion	

### Accessories Nano Spray Dryer B-90 HP




#### Nebulizers Set

Nebulizers in the sizes small, medium and large allow the finding of the perfect balance between small particles and high throughput (up to 200 mL / h).




#### Dehumidifier S-396

An efficient way to obtain constant parameters by conditioning the inlet air allowing to work with organic solvents and water mixtures in combination with the Inert Loop.



#### Inert Loop S-395


Used to work in closed mode, under nitrogen atmosphere with the BUCHI Spray Dryers. Furthermore it condenses the organic solvents and recirculates the nitrogen.



#### Aspirator


If no compressed air is available, the BUCHI Aspirator unit with inlet filter can be used to establish the required drying air flow rate. In “closed loop” mode, with the Nano Spray Dryer advanced, the Aspirator is required to build up the gas stream.

### Accessories Mini Spray Dryer S-300



#### Inert Loop S-395

Used to work in closed mode, under nitrogen atmosphere with the BUCHI Spray Dryers. Furthermore it condenses the organic solvents and recirculates the nitrogen.




#### Ultrasonic Package

Allows the Mini Spray Dryer to produce particles in the size range from 10 – 60 µm.



#### Dehumidifier

An efficient way to obtain constant parameters by conditioning the inlet air allowing you to work with organic solvents and water mixtures in combination with the Inert Loop.




#### High Performance Cyclone

Specially optimized to collect small particles in high yields from the Mini Spray Dryer.



#### Two-fluid Nozzle

The two fluid nozzle equipped with an effective nozzle cleaning mechanism and a ruby stone to guarantee reproducibility offers a high degree of flexibility.



#### Outlet Filter

Collects the residual particles from the cyclone and protects the user, the environment and the instrument. Available with a polyester deep filter and a PTFE filter membrane. Recommended for all Mini Spray Dryers.

# Prilling by Vibration

## Product Details [↗](#)




### System Portfolio & Technical Features

Encapsulator B-390 / 395 Pro [↗](#)

Particle Size	150 - 2000 $\mu\text{m}$
Particle Nature	wet (hydrogel beads / capsules)
Particle Size Distribution	uniform
Max. Sample Throughput	600 mL / h
Min. Sample Amount	5 mL
Yield	up to 100 %
Sample Viscosity	up to 300 cps
Sample Composition	aqueous organic solutions, suspensions or emulsions wax / melt
Application	agglomeration matrix encapsulation liquid core encapsulation

### Accessories




#### Dry Air Supply Cpl.

Brings dry and clean air to work with pressure bottle on the Encapsulator.



#### Alginate Powder

The alginate is tested for microencapsulation procedures and will make your lab work more reproducible.



#### Concentric Nozzle

Used for the core-shell capsule production. Includes a pulsation chamber plus a set of 7 external nozzles with high precision opening of 0.2, 0.3, 0.4, 0.5, 0.6, 0.7 and 0.9 mm.

# Melting Point


Product Details [↗](#)




## System Portfolio & Technical Features

Melting Point Model	Compliant with Pharmacopeia Methods	Automatic Detection	Sample Loader
M-560 <a href="#">↗</a>	yes	no	no
M-565 <a href="#">↗</a>	yes	yes	no
M-565 + Sample Loader M-569 <a href="#">↗</a>	yes	yes	yes


## Accessories




**Verification Kit**  
Kit of three BUCHI certified standards for verification of the Melting Point M-560 and M-565.




**Calibration Kit**  
Kit of four BUCHI certified standards for calibration of the Melting Point M-560 and M-565.




**MeltingPoint Monitor Software with License**  
Software CD, Tutorial and single PC license. For installation under Windows 7 Professional / Enterprise / Ultimate (32-bit or 64-bit, SP1), Windows 8.1 Professional / Enterprise (64-bit), Windows 10 Professional / Enterprise (64-bit).



**Sample Loader**  
Instrument for fast and efficient loading of samples into melting point capillaries.




**Printer and Keyboard**  
For documenting calibration results and melting and boiling point determinations and convenient straightforward input of parameters.




**IQOQ Documentation**  
Installation / Operational Qualifications documentation.

## Consumables



**Melting Point Capillaries**  
The precision glass capillaries ensure highly reproducible melting point determinations.




**Boiling Point Tube A Boiling Point Capillary B**  
The precision capillaries B generate perfect gas bubbles within tube A for reproducible boiling point determinations.

# NIR Spectroscopy

## Product Details - ProxiMate™ [↗](#)




### Accessories




#### Sample Cups

Glass petri dish, High Performance Sample Cup with certified glass bottom, Robust Cup with hardened glass bottom and silicone ring.



#### Large Sample Cup


For the analyses of inhomogeneous samples in Down-View Mode.



#### Transflectance Cover

For the measurement of liquid samples in Up-View mode.

### Consumables



#### PS Petri Dish

For the analyses of samples in Down-View mode.

### System Portfolio & Technical Features [↗](#)

ProxiMate™ Mode	Sample Types	Measurement Modes	Sample Cups	Optional Vis Detector <a href="#">↗</a>
Up-View	Powders, granulates, liquids, gels	Reflexion, Transflectance	Glass petri dish, High Performance Sample Cup, Robust Cup	Yes
Down-View	Powders, granulates, very viscous and non-translucent gels	Reflexion	PS petri dish, large sample cup	Yes
Dual-View	Powders, granulates, liquids, gels	Reflexion, Transflectance	All sample cups	Yes

# NIR Spectroscopy


## Product Details - NIRFlex N-500 [↗](#)




### System Portfolio & Technical Features [↗](#)

NIRFlex N-500 Measurement Cell	Sample Types	Measurement Mode
Solids	Powders, granulates, liquids, gels	Reflexion, Transflectance
Liquids	Liquids	Transmission
Fibre Optics Solids/Liquids	Powders, granulates, liquids, gels	Reflexion (only Fibre Optics Solids), Transflectance
Solids Transmittance	Tablets and capsules	Transmission


### Measurement Cells




**Solids**  
Different add-ons for different sample scenarios are available, see below.




**Liquids**  
For the temperature-controlled measurement of liquids in cuvettes or disposable vials.



**Fibre Optics Solids**  
For the direct NIR measurement of raw materials in containers




**Solids Transmittance**  
For the measurement of content uniformity in transmission of whole tablets and capsules.




**IQ/ OQ Documentation**  
Installation / Operational Qualifications documentation.


### Add-Ons for Solids Cell




**Petri Dish Add-On**  
To be used with glass petri dish, PS petri dish, the High Performance Sample Cup, or the Robust Cup (see page 31).



**XL Add-On With Iris Aperture**  
For sample cups smaller than petri dish.




**XL Add-On**  
For the measurements of samples in bags.




**Vial Add-On**  
For the consecutive measurement of 6 vials.

# NIR Spectroscopy


## Product Details - NIR-Online Process Analyzer



### Measurement Cells




**Weld-In Flange**  
For direct mount to the production line. No physical contact of the sensor with the product.



**X-Square**  
For all kind of solids and powders. Designed for continuous product flows.




**X-Cell DN50, PL1, 10bar**  
For liquid and paste-like products.




**XL-Feeder**  
For all kind of solids and powders. Designed for discontinuous product flows.


### Software



**SX-Suite**  
Intuitive user interface provides the interpretation of the instrument data and structures it to monitor, control and document the process.



**AutoCal**  
AutoCal is the most convenient tool available on the market to create and maintain calibrations.



**SX-Plus**  
Toolbox to create, update and optimize calibrations.

### System Portfolio & Technical Features

	X-One 	X-Sential™ 	Multipoint System 	PA2 
Sample Type	Solids, liquids, viscous, pasty	Solids, liquids, viscous, pasty	Solids, liquids, viscous, pasty	Solids, liquids, viscous, pasty
ATEX Certification	Dust- & Gas-Ex	General purpose, only	Dust-Ex	Dust- & Gas-Ex
When to use	<ul style="list-style-type: none"> <li>· for demanding process applications</li> <li>· for fast moving goods e.g., conveyor belt installations</li> <li>· for NIR or VIS measurements</li> <li>· for visible detections of e.g., foreign particles</li> </ul>	<ul style="list-style-type: none"> <li>· for monitoring basic parameters such as moisture, fat, or protein</li> <li>· developed and designed as the most cost-efficient choice</li> </ul>	for covering the production with up to 10 measurement points with only one spectrometer	for demanding process applications requesting an extended spectrometer range

# Core messages to our customers

## BUCHI creates added value

“Quality in your hands” is the guiding principle that shapes our philosophy and our actions. It challenges us to provide outstanding services that are precisely tailored to your needs. This means that we must stay in close contact with our customers. That is why we keep in touch and continue to work very hard to understand you and your business even better.

We help you by providing high-quality products, systems, solutions, applications and services that offer you added value. This allows you to focus entirely on your processes and your work.



### Competent

We have the technological expertise and decades of experience needed to provide competent support and work with you to continually improve our services.



### Reliable

We guarantee the quality and functionality of our equipment and will continue to help you quickly and efficiently whenever something does not operate to your satisfaction.



### Safe

By collaborating closely with you, we do everything in our power to make our products, systems, solutions, applications and services as safe as possible for people and the environment.



### Cost-effective

We strive to create a high level of economic benefit and maximum added value for you.



### Global

As an international family-owned business with own subsidiaries and qualified distributors, we have a presence wherever you are located.



### Easy

We support you by providing carefully designed solutions as well as instruments and systems that are easy to operate.



### Sustainable

We support environmentally friendly processes and manufacture products that have a long service life. We utilize advanced technologies to leave the smallest environmental footprint possible.

---

We are represented by more than 100 distribution partners worldwide.  
Find your local representative at:

[www.buchi.com](http://www.buchi.com)

Quality in your hands

