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







	注文番号	イラスト
Heating Bath B-301, 110V, not incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95° C. Meant to be used with water as heating media.	11B301002	
Heating Bath B-301, 230V, not incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95° C. Meant to be used with water as heating media.	11B301001	
Heating Bath B-301, 110V, incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95° C. Meant to be used with water as heating media.	11B301102	

<emphasis type="Accessories" />

	注文番号	イラスト
Heating Bath B-301, 230V, incl. Base B-300 For up to 1 Liter evaporating flasks. For temperatures up to 95 ° C. Meant to be used with water as heating media.	11B301101	
Heating Bath B-305, 110V, not incl. Base B-300, 110V For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media.	11B305002	
Heating Bath B-305, 230V, not incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media.	11B305001	
Heating Bath B-305, 110V, incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media.	11B305102	
Heating Bath B-305, 230V, incl. Base B-300 For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media.	11B305101	
Heating Bath B-305, 110V, for R-21x, incl. Base For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media. Included: Adapter for Rotavapor® R-21x and R-II.	11B305112	
Heating Bath B-305, 230V, for R-21x, incl. Base For up to 5 Liter evaporating flasks. For temperatures up to 220 ° C. Meant to be used with water or oil as heating media. Included: Adapter for Rotavapor® R-21x and R-II.	11B305111	

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	注文番号	イラスト
Dewar accessory, set. Incl. Dewar container and adapter For sample preparation in freeze drying. To be used with dry ice and ethanol / isopropanol / acetone or separately with liquid nitrogen. Compatible with Heating Bath B-300 (not plugged in during operation) and Rotavapor® R-300.	11066662	
Dewar container, not incl. adapter For sample preparation in freeze drying. To be used with dry ice and ethanol / isopropanol / acetone or separately with liquid nitrogen. To be used only with adapter, Heating Bath Base B-300 (not plugged-in during operation) and Rotavapor R-300.	11066645	
Lid. For Heating Bath B-305 To save energy and to minimize loss of water when heating bath is not in use.	11059500	
Protection shield, Operator protection for heating bath For Heating Bath B-301 and B-305. Material: 4 mm thick PMMA.	11061402	
Splash shield, Operator protection against splashes For Heating Bath B-301 and B-305. Material: 6 mm thick polycarbonate.	11064768	
Protective cover. For Heating Bath B-300 base Used for protection of the heating bath display.	11064111	

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


注文番号
Heating bath balls. 450 pcs., PP, Ø10 mm
To reduce energy consumption of heating bath and for less evaporation of the heating medium. For temperatures up to 100 ° C.

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








注文番号
Tubing. FEP, Ø6/8 mm, transparent, per m
Use: Vacuum, cooling media.

注文番号	
Tubing. Natural rubber, Ø6/16 mm, red, per m	017622
Use: Vacuum.	
Tubing. Nylflex, PVC-P, Ø8/14 mm, transparent, per m	004113
Tubing. PTFE, Ø4.7/5.5 mm, transparent, 330 mm	000646
Use: To introduce solvent into evaporating flask during distillation.	
Tubing. PTFE, Ø4.7/5.5 mm, transparent, 460 mm	000643
Use: To introduce solvent into evaporating flask during distillation.	
Tubing. Silicone, Ø6/9 mm, transparent, per m	004133
Use: Cooling media.	
Tubing. Synthetic rubber, Ø6/13 mm, black, per m	11063244
Use: Vacuum.	

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

注文番号	イラスト
StatusLight cpl., incl. communication cable	11068959
Indicates the status of the instrument (instrument is ready to use, has an error or is in operation).	
BUCHI Bluetooth® ドングル、本機とスマートフォンをBluetooth®で接続します。	11067770
	
AutoDest sensor. Incl. cap nut, seal GL14	11059225
For automatic distillation. Measures temperature of cooling media and the vapor temperature. Vacuum is adjusted according to cooling capacity of condenser. Meant to be used with the Interface I-300 / I-300 Pro and glass assembly V, HP or S.	
Condensate trap.	11062955
Condenser V and HP, grey, TPE, hose connector Ø8 mm. Collects and drains condensate which may accumulate on the condenser.	
Condensate trap.	11061985
Condenser V and HP, green, TPE, hose connector Ø8 mm. Collects and drains condensate which may accumulate on the condenser.	
Cooling water valve. 24VAC	031356
Valve opens cooling water feed during distillation. Meant to be used with a vacuum controller/interface.	

注文番号	イラスト
円筒ろ紙セット25個、43 x 118 mm、セルロース ソックスレー抽出装置用。 作業容量：150 mL	018106 
Flask holder. EPDM, slip free Holder for round-bottom flasks (50-5000 mL).	048618 
Flask holders, set. 5 pcs., EPDM, slip free Holder for round-bottom flasks (50-5000 mL).	11059916 
Foam sensor. Incl. holder Prevents sample from foaming into the condenser by aeration of system temporarily. Meant to be used with the Interface I-300 / I-300 Pro.	11061167 
IQ/OQ. R-300 en Official BUCHI documents.	11063608 
Back-feeding valve, cpl.	048100
Lab grease. Glisseal 40 (30 g) To grease joints in order to increase tightness of system.	048197 
Level sensor. Incl. O-ring, cable, rubber band Prevents an overflow in the receiving flask of condenser/secondary condenser. Placed on receiving flask. Meant to be used with the Interface I-300 / I-300 Pro.	11060954 
Repeating OQ. R-300 en Official BUCHI documents.	11063609 
Power plug fixation. For R-300 Used to fixate the power plug into Rotavapor®.	11063710 
Set earthquake fastener. Lug for backside To fasten device on lab bench.	11062386 

	注文番号	イラスト
Thread adapter. For vacuum connection, SVL22/ GL14	040871	
Use: To switch vacuum connection from a Rotavapor® to Syncore® vacuum distributor (037686).		
Vacuum connection. Incl. drain valve, hose barb Ø 8 mm, BJ35/20	001006	
Connection piece for aeration of the system, placed between condenser and receiving flask.		
Vacuum valve. Magnetic valve, 24V/4W, Mini-DIN, 1.5 m	11060706	
Flow valve without flask, meant to be used with a centralized vacuum source or an unregulated vacuum pump. Meant to be used with the Interface I-300 / I -300 Pro.		
Vapor temperature sensor. Incl. cap nut, seal GL14	11060707	
Measures the vapor temperature inside the system. Meant to be used with the Interface I-300 / I-300 Pro.		
Water jet pump. Plastic	002913	
Used when tap water is used to generate vacuum.		
Water regulation nozzle. Flow regulator, incl. hose clamp, sieve	011606	
Used when tap water is used to generate vacuum. Reduces water consumption.		
Woulff bottle. 125 mL, P+G, incl. holder	11059905	
For trapping particles and droplets and for pressure equalization.		
Woulff bottle. 3-neck, 800 mL, P+G.	025519	
For trapping particles and droplets and for pressure equalization.		
Manometer with needle valve. Incl. vacuum gauge, not incl. holder	047391	

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

	注文番号	イラスト
Glass assembly A. Diagonal condenser, 1500 cm ² For standard distillations, used with recirculating chiller or tap water. Features: Cleaning hole on top. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048168	
Glass assembly A. Diagonal condenser, 1500 cm ² , P+G For standard distillations, used with recirculating chiller or tap water. Features: Cleaning hole on top. Content: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048169	
Glass assembly BF. Vertical condenser, 1500 cm ² For simultaneous back-feeding and foam prevention, used with recirculating chiller or tap water. Features: Large cleaning hole on top (SVL22) and on the side (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing, 048100 valve. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder. Use of foam sensor and back-feeding possible only with 29/32 vapor duct.	11074661	

	注文番号	イラスト
Glass assembly BF. Vertical condenser, 1500 cm ² P+G For simultaneous back-feeding and foam prevention, used with recirculating chiller or tap water. Features: Large cleaning hole on top (SVL22) and on the side (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing, 048100 valve. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder. Use of foam sensor and back-feeding possible only with 29/32 vapor duct.	11074662	
Glass assembly BY. Vertical intensive condenser, reflux, 1300 cm ² , P+G For particularly efficient condensation with double jacket condenser as well as for reflux distillations with distributor and shut-off valve. Used with recirculating chiller or tap water. Features: Double jacket condenser. Additional joint on top of the condenser for flexible expansion (SJ29/32). Connection for vapor temperature sensor. Content: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing, distributor. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048297	
Glass assembly C. Cold trap, 500 cm ² For distillations of low-boiling point solvents. No need of cooling liquid, but e.g. dry ice. Content: 1 L receiving flask, ball joint clamp, stopcock, feeding tubing. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	040640	



注文番号	イラスト
Glass assembly C. Cold trap, 500 cm ² , P+G-LT For distillations of low-boiling point solvents. No need of cooling liquid, but e.g. dry ice. Application temperature: -70 to 40 ° C. Content: 1 L receiving flask (P+G-LT), ball joint clamp, stopcock, feeding tubing. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	040642
Glass assembly CR. Cold trap reflux, 500 cm ² For reflux distillations with distributor and shut-off valve as well as for distillation of low-boiling point solvents. No need of cooling liquid, but e.g. dry ice. Content: 1 L receiving flask, ball joint clamp, stopcock, feeding tubing, distributor. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048292
Glass assembly CR. Cold trap reflux condenser, 500 cm ² , P+G-LT For reflux distillations with distributor and shut-off valve as well as for distillation of low-boiling point solvents. No need of cooling liquid, but e.g. dry ice. Application temperature: -70 to 40 ° C. Content: 1 L receiving flask, ball joint clamp, stopcock, feeding tubing, distributor. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048293
Glass assembly E. Expansion condenser, 1500 cm ² Optimised for distillations of foaming substances and bumping sample solutions. Used with recirculating chiller or tap water. Features: Descending condenser with expansion vessel. Cleaning hole on top. Connection for vapor temperature sensor. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing, holder. Not inclusive: Evaporating flask, vapor duct, vacuum gasket.	11061112

	注文番号	イラスト
Glass assembly E. Expansion condenser, 1500 cm ² , P+G Optimised for distillations of foaming substances and bumping sample solutions. Used with recirculating chiller or tap water. Features: Descending condenser with expansion vessel. Cleaning hole on top. Connection for vapor temperature sensor. Content: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing, holder. Not inclusive: Evaporating flask, vapor duct, vacuum gasket.	11061113	
Glass assembly S. Vertical condenser reflux, 1500 cm ² For reflux distillations with distributor and shut-off valve. Used with recirculating chiller or tap water. Features: Cleaning hole on top. Connection for vapor temperature sensor. Inner glass tube. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing, distributor. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048290	
Glass assembly S. Vertical condenser reflux, 1500 cm ² , P+G For reflux distillations with distributor and shut-off valve. Used with recirculating chiller or tap water. Features: Cleaning hole on top. Connection for vapor temperature sensor. Inner glass tube. Content: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing, distributor. Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.	048291	


注文番号	イラスト
Glass assembly V. Vertical condenser, 1500 cm ²	11062432
For standard distillations, used with recirculating chiller or tap water.	
Features: Large cleaning hole on top (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Content: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing.	Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.
Glass assembly V. Vertical condenser, 1500 cm ² , P+G	11062433
For standard distillations, used with recirculating chiller or tap water.	
Features: Large cleaning hole on top (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Content: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing.	Not inclusive: Evaporating flask, vapor duct, vacuum gasket, holder.
Glass assembly HP, Incl. condenser, supplementary set and holder	11066561
High performance condenser, 3000 cm ² , for faster process, greater distillation rates and fewer emissions of solvents, used with recirculating chiller or tap water. Features: Large cleaning hole on top (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Contents: 1 L receiving flask, required tubings, ball joint clamp, stopcock, feeding tubing, holder.	
Not included: Evaporating flask, vapor duct, vacuum gasket. Compatible with Rotavapor® R-300, R-215 and R-210.	


	注文番号	イラスト
<p>Glass assembly HP, Incl. condenser P+G, supplementary set and holder</p> <p>High performance condenser, 3000 cm², P+G, for faster process, greater distillation rates and fewer emissions of solvents, used with recirculating chiller or tap water. Features: Large cleaning hole on top (SVL22). Connection for vapor temperature sensor. Inner glass tube. Inner glass plate. Contents: 1 L receiving flask (P+G), required tubings, ball joint clamp, stopcock, feeding tubing, holder.</p> <p>Not included: Evaporating flask, vapor duct, vacuum gasket. Compatible with Rotavapor® R-300, R-215 and R-210.</p>	11066562	

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
	注文番号	イラスト
<p>Condenser holder. For glass assembly V/HP/C/S/ CR/BY</p> <p>Content: Holding rod, rubber band, cross sleeve.</p>	048180	
<p>Condenser holder. For Glass assembly E</p> <p>Content: Holding rod, bride, cross sleeve.</p>	048181	
<p>Extraction unit Soxhlet, 200 mL, incl. extraction part, reduction part</p> <p>For Soxhlet extraction applications with a Rotavapor® R-300. Meant to be used with glass assembly S. Condenser holder is included.</p>	011744	
<p>Extraction unit Soxhlet, 500mL, incl. extraction part, reduction part</p> <p>For Soxhlet extraction applications with a Rotavapor® R-300. Meant to be used with glass assembly S. Condenser holder is included.</p>	011745	

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






	注文番号
Glass, SJ24/40, 500 mL With 4 indents for better mixing/drying.	011579
Glass, SJ29/32, 500 mL With 4 indents for better mixing/drying.	000452
Glass, SJ24/40, 1000 mL With 4 indents for better mixing/drying.	000420
Glass, SJ29/32, 1000 mL With 4 indents for better mixing/drying.	000453
Glass, SJ24/40, 2000 mL With 4 indents for better mixing/drying.	011580
Glass, SJ29/32, 2000 mL With 4 indents for better mixing/drying.	000454

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	注文番号	イラスト
Beaker flask, Drying, bayonet type, flat bottom, SJ24/40, 500 mL Bayonet mount type Ø75 mm. With 4 indents for better mixing/drying. Working volume 150 mL. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.	11063159	
Beaker flask, Drying, bayonet type, flat bottom, SJ24/40, 1500 mL Bayonet mount type Ø110 mm. With 4 indents for better mixing/drying. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.	11063161	
Beaker flask, Drying, bayonet type, flat bottom, SJ29/32, 500 mL Bayonet mount type Ø75 mm. With 4 indents for better mixing/drying. Working volume 150 mL. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.	11063158	

	注文番号	イラスト
Beaker flask, Drying, bayonet type, flat bottom, SJ29/32, 1500 mL	11063160	
Bayonet mount type Ø110 mm. With 4 indents for better mixing/drying. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, flat bottom, SJ24/40, 500 mL	11063155	
Bayonet mount type Ø75 mm. Working volume 150 mL. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, flat bottom, SJ24/40 1500 mL	11063157	
Bayonet mount type Ø110 mm. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, flat bottom, SJ29/32, 500 mL	11063154	
Bayonet mount type Ø75 mm. Working volume 150 mL. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, flat bottom, SJ29/32 1500 mL	11063156	
Bayonet mount type Ø110 mm. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, round bottom, SJ24/40, 1500 mL	11065719	
Bayonet mount type Ø110 mm. Use: Easy cleaning and collection of solid residue. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, round bottom, SJ29/32, 1500 mL	11065718	
Bayonet mount type Ø110 mm. Use: Easy cleaning and collection of solid residue. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		
Beaker flask, Evap., bayonet type, pear shape, SJ24/40, 4000 mL	11065691	
Bayonet mount type Ø110 mm. Appropriate for natural product applications. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		

	注文番号	イラスト
Beaker flask, Evap., bayonet type, pear shape, SJ29/32, 4000 mL	11065690	
Bayonet mount type Ø110 mm. Appropriate for natural product applications. Content: 1 coupling piece, 1 seal, 1 flask and 1 beaker fastener.		

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
	注文番号	イラスト
Receiving flask with drain valve. Glas, BJ35/20, 1000 mL, P+G, PTFE	036919	

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	注文番号
Glass, BJ35/20, 50 mL	000421
Glass, BJ35/20, 100 mL	000422
Glass, SJ35/20, 250 mL	000423
Glass, SJ35/20, 250 mL, P+G	11060907
Glass, BJ35/20, 250 mL, P+G-LT	11060908
Application temperature: -70 to 40 ° C.	
Glass, SJ35/20, 500 mL	000424
Glass, SJ35/20, 500 mL, P+G	025264
Glass, BJ35/20, 500 mL, P+G-LT	040774
Application temperature: -70 to 40 ° C.	
Glass, BJ35/20, 1000 mL, P+G	020728
Glass, BJ35/20, 1000 mL, P+G-LT	040775
Application temperature: -70 to 40 ° C.	
Glass, BJ35/20, 2000 mL	000426
Glass, BJ35/20, 2000 mL, P+G	025265
Glass, BJ35/20, 3000 mL, P+G-LT	040777
Application temperature: -70 to 40 ° C.	

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
	注文番号	イラスト
For condenser C/CR, glass, SJ18.8/38	040628	
For aeration of the system. For cold trap outer part.		
Standard, glass, SJ18.8/38	040627	
For aeration of the system.		
Stopcock, Analytic PTFE/25% glass fiber, SJ18.8/38	11069607	
For feeding of solvents and aerating the system. Less cross-contamination compared to standardstopcock. For applications where grease should be avoided.		
Content: PTFE stopcock (no tubing included).		

	注文番号	イラスト
PTFE, incl. 3-way valve For feeding of solvents and aeration of the system. For applications when grease should be avoided. Used instead of standard-stopcock (040627). Content: Inlet tubing 300mm, backfeed tubing 600mm, cap nut GL10.	11058814	

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	注文番号
Glass, SJ24/40, incl. 100 mL flask (5pcs)	011575
Glass, SJ24/40, incl. 50 mL flask (5pcs)	011574
Glass, SJ29/32, incl. 100 mL flask (5pcs)	001333
Glass, SJ29/32, incl. 20 mL cyl. flask (12pcs)	001335
Glass, SJ29/32, incl. 20 mL cyl. flask (20pcs)	001336
Glass, SJ29/32, incl. 20 mL cyl. flask (6pcs)	001334
Glass, SJ29/32, incl. 50 mL flask (5pcs)	001332
Evaporating flask, For distillation spider, cylindric, SJ14/23, 20 mL	000477

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	注文番号	イラスト
Intermediate piece with valve, Incl. drain valve, hose connection Ø8 mm, BJ35/20	11063430	
Connection piece with 3-way valve, for aeration of the system. Placed between condenser and receiving flask. Allows to remove and empty the receiving flask during the evaporating process.		
Set Rotavapor® connection Combining the Multivapor with the Rotavapor® requires a T-shaped glass connector for the condenser of the rotary evaporator. The T-piece consists of two spherical joints for the condenser assembly and a SVL 22 joint for the vacuum tube. The length of the tube is 400 mm.	048740	

注文番号	イラスト
Vacuum connection. Incl. drain valve, hose barb Ø 8 mm, BJ35/20	001006
Connection piece for aeration of the system, placed between condenser and receiving flask.	

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

注文番号
Vacuum seal VS26, PTFE base, NBR O-ring, FDA-complaint

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注文番号
For cap nut GL14 to FEP, EPDM
Set. 10pcs, for hose barbs GL14, EPDM, black
Set. 10pcs, for hose barbs GL14, FPM, green
Set. 20pcs, for hose barbs GL14, silicone, red

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注文番号	イラスト
Hose barb. Bent, GL14, incl. silicone seal	018916
Hose barbs, set. 2 pcs, bent (1), straight (1), GL14, silicone seal Content: Hose barbs, cap nuts, seals	041939
Hose barbs. set. 3 pcs, bent, GL14, silicone seal Content: Hose barbs, seals.	041987
Hose barbs, set. 4 pcs, bent GL14, silicone seal Content: Hose barbs, cap nuts, seals	037287
Hose barbs. set. 4 pcs, bent, GL14, EPDM seal Content: Hose barbs, cap nuts, seals.	043129
Hose barbs, set 4 pcs, bent, GL14, FPM seal Content: Hose barbs, cap nuts, seals.	040295

	注文番号	イラスト
Hose barbs. set. 4 pcs, straight, GL14, EPDM seal Content: Hose barbs, cap nuts, seals.	043128	
Hose barbs. set. 4 pcs, straight, GL14, FPM seal Content: Hose barbs, cap nuts, seals.	040296	
Hose barbs, set. 4 pcs, straight, GL14, silicone seal Content: Hose barbs, cap nuts, seals	037642	
Hose barbs. set. 6 pcs, bent (4), straight (2), GL14, silicone seal Content: Hose barbs, cap nuts, seals.	038000	
その他の消耗品		
Cap nut. Screw cap with hole GL10	023875	
Cap nuts, set. 10 pcs, screw cap with hole GL14, seal EPDM Content: Cap nuts, seals	041999	
Cap nuts. set. 10 pcs, screw cap with hole, GL14	041956	
Hose, set. Incl. GL14 set, FEP tube 11061756 Used to connect VacuBox and Valve unit/Woulff bottle when both are fixed on the Rotavapor® R-300.	11065373	
O-ring. FKM, 40.9/2.6 mm, for reflux insert	048078	
O-ring. FPM, for cap nut GL10, Ø3.00/2.70 mm	023900	







	注文番号	イラスト
Screw caps. set. 5 pcs, GL14	040624	
Screwing set, Hose barbs and seals GL14 Content: Hose barbs bent 4 pcs, hose barbs straight 1 pc, cap nuts 4 pcs, screw caps 2 pcs, seals EPDM 3 pcs.	11061921	

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
	注文番号	イラスト
Beaker flask, Drying, screw joint, flat bottom, glass, 500 mL With 4 indents for better mixing/drying. Compatible with a beaker fastener (034771) Ø75 mm.	034770	
Beaker flask, Drying, screw joint, flat bottom, glass, 1500 mL With 4 indents for better mixing/drying. Compatible with a beaker fastener (034139) Ø110 mm.	034267	
Beaker flask, Evap., screw joint, flat bottom, glass, 500 mL Compatible with a beaker fastener (034771) Ø75 mm.	034769	
Beaker flask, Evap., screw joint, flat bottom, glass, 1500 mL Compatible with a beaker fastener (034139) Ø110 mm.	034140	
Beaker flask, Drying, bayonet type, flat bottom, glass, 500 mL With 4 indents for better mixing/drying. Working volume 150 mL. Compatible with a beaker fastener (11059810) Ø75 mm.	11059268	
Beaker flask, Drying, bayonet mount, flat bottom, glass, 1500 mL With 4 indents for better mixing/drying. Compatible with a beaker fastener (11059851) Ø110 mm.	11059269	

	注文番号	イラスト
Beaker flask, Evap., bayonet type, flat bottom, glass, 500 mL Working volume 150 mL. Compatible with a beaker fastener (11059810) Ø75 mm.	11059185	
Beaker flask, Evap., bayonet type, round bottom, glass, 1500 mL Compatible with a beaker fastener (11059851) Ø110 mm.	11065716	
Beaker flask, Evap., bayonet type, flat bottom, glass, 1500 mL Compatible with a beaker fastener (11059851) Ø110 mm.	11059186	
Beaker flask, Evap., bayonet type, pear shape, glass, 4000 mL Compatible with a beaker fastener (11059851) Ø110 mm.	11065689	
Coupling piece, Bayonet mount type, glass, SJ 24/40, Ø75 mm For 500 mL beaker flasks.	11059362	
Coupling piece, Bayonet mount type, glass, SJ29/32, Ø75 mm For 500 mL beaker flasks.	11059043	
Coupling piece, Bayonet mount type, glass, SJ 24/40, Ø110 mm For 1500 mL round and flat bottom and for 4000 mL pear shape beaker flasks.	11059363	
Coupling piece, Bayonet mount type, glass, SJ29/32, Ø110 mm For 1500 mL round and flat bottom and for 4000 mL pear shape beaker flasks.	11059187	
Distillation spider. Glass, SJ29/32 (for 20 flasks), not incl. flasks	000476	
Distillation spider. Glass, SJ24/40 (for 5 flasks), not incl. flasks	000466	
Distillation spider. Glass, SJ29/32 (for 12 flasks), not incl. flasks	000475	



	注文番号	イラスト
Distillation spider. Glass, SJ29/32 (for 5 flasks), not incl. flasks	000471	
Distillation spider. Glass, SJ29/32 (for 6 flasks), not incl. flasks	000474	
Distribution head. For assembly S/CR/BY, P+G, incl. cap nut GL10	040658	
Distribution head. For glass assembly S/CR/BY, incl. cap nut GL10	040657	
Distribution head. For glass assembly E, incl. cap nut GL10	040661	
Distribution head. For glass assembly E, P+G, incl. cap nut GL10	040663	
Expansion vessel. For Glass assembly E, glass, incl. screw cap GL14	001002	
Expansion vessel. For Apparatus E, glass, P+G, inc. GL14 union nut	033507	
U-tube. For assembly E, Glass, incl. screw cap GL14	001003	
U-tube. For assembly E, Glass, P+G incl. screw cap GL14	033508	
Cold finger. For condenser C/CR	000672	





	注文番号	イラスト
Vacuum connection. For glass assembly BY, vacuum piece	035026	
To close the top hole (NS29/32) of the condenser if no expansion is used. Openings with screw cap SVL15 (on top) and with cap nut GL14 (sideward).		

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	注文番号
Gasket, For 1.5 L and 4 L beaker flask, Ø100 mm, incl. O-ring	034152
PTFE coated gasket is used to seal the beaker flask.	
For beaker flask 500 mL, Ø68 mm, incl. O-ring	034772
PTFE coated seal for beaker flask.	
Set. For condenser C/CR, O-ring EPDM/PTFE	027462
Seal between cold trap outer part and cold finger.	
Set. For Woulff bottle/valve unit	045544
Content: O-ring (048406, 5pcs), O-ring (047169, 5pcs).	
Set. For Woulff bottle/valve unit	045545
Content: O-ring (048406, 5pcs), O-ring (047169, 5pcs), spare parts for valve unit (047258), seal (047165).	
Seal, SVL sealing for foam sensor	11058958

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








	注文番号	イラスト
Condenser A. Diagonal condenser, 1500 cm ²	046988	
Condenser A. Diagonal condenser, 1500 cm ² , P+G	047976	
Condenser BF. Vertical condenser, 1500 cm ²	11074616	







	注文番号	イラスト
Condenser BF. Vertical condenser, 1500 cm ² , P+G	11074617	
Condenser BY. Double jacket, 1300 cm ² , P+G	035040	
Cold trap mantle. For condenser C	040641	
Condenser C. Cold trap, 500 cm ²	040645	
Condenser C. Cold trap, 500 cm ² .P+G-LT Application temperature: -70 to 40 ° C.	040646	
Condenser CR. Cold trap reflux, outer part	011228	
Condenser CR. Cold trap reflux, outer part, P+G-LT Application temperature: -70 to 40 ° C.	025614	
Condenser CR. Cold trap reflux, 500 cm ²	011511	
Condenser CR. Cold trap reflux, 500 cm ² , P+G-LT Application temperature: -70 to 40 ° C.	033478	
Condenser E. Expansion, 1500 cm ²	010183	
Condenser E. Expansion, 1500 cm ² , P+G	033511	
Condenser S. Vertical condenser reflux, 1500 cm ²	040651	

	注文番号	イラスト
Condenser S. Vertical condenser reflux, 1500 cm ² , P+G	040653	
Condenser V. Vertical condenser, 1500 cm ²	11059368	
Condenser V. Vertical condenser, 1500 cm ² , P+G	11059369	
Condenser HP, High performance condenser, 3000 cm ² Compatible with Rotavapor® R-300, R-215 and R-210.	11065785	
Condenser HP, High performance condenser, 3000 cm ² , P+G Compatible with Rotavapor® R-300, R-215 and R-210.	11065786	

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	注文番号	イラスト
Ball joint clamp. For BJ35/20	003275	
To fasten receiving flask on condenser/secondary condenser.		
Beaker fastener, Connection with bayonet mount, for 1.5 and 4 L	11059851	
For 1500 mL round and flat bottom and for 4000 mL pear shape beaker flasks.		
Beaker fastener, Connection with bayonet mount, for 500 mL	11059810	
For beaker flask 500 mL.		
Capacitive switch, For level sensor	11055914	
Combi-Clip, Not incl. vapor duct	11059770	
Combi-Clip with snap lock mechanism to fasten evaporating flask on vapor duct.		

	注文番号	イラスト
Combi-Clip, set, 50pcs, not incl. vapor duct Combi-Clip with snap lock mechanism to fasten evaporating flask on vapor duct.	11063308	
Cross sleeve (1 pc.)	027344	
Deflector rings, set. Drain sleeve, EPDM, black, Ø41/55 mm Drains off condensate which may accumulate at the condenser.	040822	
Drain disc. PTFE, Ø5.1/14 mm	040625	
Flange screwed connection. For Glass assembly V/HP/C/A Content: Flange nut, pressure spring.	11062387	
Flange screwed connection-reflux. For Glass assembly E/S/BY/CR Content: Flange nut, reflux part, pressure spring, O-ring FKM.	048080	
Fuses, set. 10pcs, 6.3AT, 20 mm, Ø5 mm	022561	
Fuses, set. 10pcs, 12.5AT, 20 mm, Ø5 mm	047939	
Handle, For top cover heating bath (B-305, 11059500) Handle with top cover is used for covering the heating bath after the experiment is over.	11059502	
Heating Bath B-300 Base. For Heating Bath B-301 and B-305, 110 V	11B300102	
Heating Bath B-300 Base. For Heating Bath B-301 and B-305, 230 V	11B300101	

	注文番号	イラスト
Holder. Clamp for glass assembly E	048125	
Holder, set, For fixing the Woulff bottle onto R-300 Content: Holder Woulff bottle and 2 M4x8 screws.	11065756	
Navigation knob. Ø30mm, green/grey, incl. knob case, cover For Rotavapor® R-300, Heating Bath B-300 Base, Interface I-300/I-300 Pro and Recirculating Chiller F-305/F-308/F-314.	11059157	
Receiving vessel. For Woulff bottle, GL40, 125 mL, P+G	047233	
Rubber band. To fasten condenser on holder	032013	
Sieve. For water control valve (011606), Ø18mm	011514	
Shut-off valve. For distribution head of glass assembly S/CR/BY For reflux reactions, for interruption of condensation.	040626	
Support foam sensor, Incl. clamping nut, support ring To fasten foam sensor on glass assembly	11059024	
T-piece. Incl. cap nut (3Stk) GL14	038190	
Tubing. PTFE, Ø3/4 mm, white, 600 mm Use: Feeding.	028096	
円筒ろ紙セット25個、43 x 118 mm、セルロース ソックスレー抽出装置用。 作業容量：150 mL	018106	

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Introduction Training

Educating Rotavapor® Users

The following topics shall be trained as part of the “Introduction Training” after the installation:

1. **Theoretical basis**
 - Reviewing fundamental principles of distillation under vacuum
 - Explanation of the Rotavapor® working principle
 - Understanding the impact of the Δ20 rule, pressure adjustment, condenser loading, rotation speed and flask size
2. **Rotavapor® R-300**
 - Setting the heating bath temperature, condenser angle and rotation without the interface
 - Locking the heating bath and chiller temperature
 - Setting and deleting a depth stop
 - Proper Combi-clip use (mounting and removing the evaporating flask, removing the vapor-duct)
 - Correctly starting and stopping the distillation process
 - Adding solvents during distillation
3. **Interface I-300 / I-300 Pro**
 - Setting the pressure, rotation, lift height and heating and cooling temperatures
 - Implementing operating modes (i.e. Manual, Timer, Continuous, AutoDest, Drying and Methods)
 - Utilizing the solvent library and favorites feature
 - Remote monitoring using the Rotavapor® App
 - Adjusting the hysteresis and maximum pump output
4. **Safety features**
 - Reviewing prohibited uses and potential dangers as described in the operation manual
 - Familiarization with P+G coating, protection and splash shield, safety stop in case of power failure, temperature limit setting, overpressure prevention, heating bath overheat protection
5. **Maintenance and cleaning**
 - Cleaning the heating bath and the vacuum pump
 - Cleaning the condenser, vapor duct, Woulff bottle, receiving and evaporating flask
 - Preventive maintenance
 - Using the leak test, finding and removing a leak
6. **Convenient additional functions**
 - Automatic foam removal (foam sensor), automatic vacuum regulation/control (AutoDest)
 - Additional safety add-ons (e.g. level sensor)

Account Details

Account Name:

Account Address:

Location (e.g. Building, Department, Lab):

Optional: Account Number (CRM):

Attendees

The people listed below have attended the “Introduction Training” for Rotavapor® users.

No.	Surname	Name	Signature	Rating*
1				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
2				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
3				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
4				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
5				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
6				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
7				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
8				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
9				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤
10				<input type="radio"/> ① <input type="radio"/> ② <input type="radio"/> ③ <input type="radio"/> ④ <input type="radio"/> ⑤

*Rate your satisfaction with the provided training. (①: very unsatisfied; ⑤: very satisfied)

Trainer

The following person confirms that he is qualified to perform the “Introduction Training” (certificate is available upon request) and covered all applicable topics mentioned on the front page.

Trainer Company Name:

Full Name of Trainer:

Date of Training:

Optional: Case Number (CRM):

Signature of Trainer:

Advice to trainer:

Hand over the original document to the customer and email a picture of this page to
InstallFeedback@buchi.com. Please use the Rotavapor® R-300 serial number as the subject line.

BÜCHI Labortechnik AG

Meierseggstrasse 40

9230 Flawil 1

Switzerland

T +41 71 394 63 63

F +41 71 394 64 64

We are represented by more than 100 distribution partners worldwide.
Find your local representative at:

www.buchi.com

Quality in your hands