Guide for electrical installation

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Introduction

EN About this guide

This guide contains information regarding the electrical installation of single-phase and three-phase devices. It is addressed to electricians or persons with similar expert knowledge.

Installers must read and understand the guide before installation. The installer should conform to all safety precautions in the document and local regulations when installing a Buchi device. Information about the warning message concept can be found in the operation manual, section "Safety".

No technical modifications may be made to the device without the prior written agreement of Buchi. Unauthorized modifications may affect the system safety and result in accidents.

This document is copyright. Information from it may not be reproduced, distributed or used for competitive purposes, nor made available to third parties.

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EN General safety rules

Quality and safety of an electrical installation

Quality and safety will be assured only if the following control procedures are respected:

- The electrical installation is conform to all applicable standards and regulations
- The electrical installation equipment comply with applicable standards and regulations
- Periodic checks of the installation /the installed device are performed according to applicable standards and regulations

Basic installation requirements and tests

Prior to installation some general electrical tests and visual inspections of the power supply network must be carried out.

These tests are made according to local regulations, which may differ slightly from one country to another. However, the basic principles are common and aim at the observance of basic safety rules in the design and realization of the installation.

A standard set of tests includes:

- Insulation tests of all cable and wiring conductors of the supply installation, between phases and between phases and ground
- Continuity and conductivity tests of protective, equipotential and ground-bonding conductors
- Verification of allowable load per circuit
- Where appropriate, verification that all exposed- and extraneous metallic parts such as distributor boxes are properly grounded
- Clearance check according to the installation requirement data of the device manual

These tests and checks are basic but not exhaustive for the majority of installation sites! Apply to local regulations and standards for more information.

 	 Danger
	Risk of death or fatal injury by electric current. • Do not install the device when missing or ineffective protective features are detected • Only install and operate the device at power networks and in environments that agree with the specifications given on the type plate • Disconnect power supply before removing cover • Cover may only be opened by instructed persons

NOTE

IEC 60364-6 and related standards are based on an international consensus for such tests, intended to cover all the safety measures and approved installation practices required in most cases.

Installation

Electrical supply

The used mains circuit has to:

- Provide the voltage and frequency that is given on the type plate of the device
- Be equipped with adequate fusage and electrical safety measures, in particular proper grounding

See also technical data of all components regarding the different system requirements and maximum power consumption!

NOTE

Switch disconnectors must meet IEC 60947-1 and IEC 60947-3 requirements. Such devices must be clearly labeled and accessible at any time.

 	Danger Risk of death or fatal injury by electric current. • Supply wires must not be under voltage prior to and during installation • Do not install the device when missing or ineffective protective features are detected • Only install and operate the device at power networks and in environments that agree with the specifications given on the type plate
 	Notice Risk of device damage by wrong mains supply. • External mains supply must meet the voltage given on the type plate • Check for proper grounding

Stationary installation – without power plug

To cut off the power, stationary installations must be equipped with an external power switch. This switch must be installed in the immediate vicinity of the installed device for fast access.

The switch must be accessible at any time during operation! The switch must be labelled prominently as a power switch!

At installation:

- Install a strain relief at the end of flexible cables
- Flexible wire ends must be equipped with a wire end ferrule

NOTE

- Additional electrical safety measures such as residual current breakers may be necessary to meet local laws and regulations!
- The grounding wire must be at least 30 mm longer, compared to all other live wires to ensure safe grounding in case of a wire fixture or strain relief failure. See drawings!

Mobile installation – equipped with power plug

Mobile installations are equipped with a power plug.

Single-phase devices

To cut off the power for single-phase devices, the socket must be installed in the immediate vicinity (max. 3 m supply cable length) of the device for fast access. The plug must not be secured with any kind of fixture!

The plug must be unpluggable at any time during operation!

Multi-phase devices

To cut off the power for multi-phase devices, an external power switch must be installed for the supply socket. This switch must be placed in the immediate vicinity of the installed device for fast access.

The switch must be accessible at any time during operation! The switch must be labelled prominently as a power switch!

At installation:

- Install a power plug with strain relief
- Use high-quality plugs (e.g. solid metal conductors only)
- Flexible wire ends must be equipped with a wire end ferrule

NOTE


- *Additional electrical safety measures such as residual current breakers may be necessary to meet local laws and regulations!*
- *External connections and extension lines must be provided with a grounded conductor lead (3-pole couplings, cord or plug equipment). All used power cords must meet the input power requirements.*
- *The grounding wire must be at least 30 mm longer, compared to all other live wires to ensure safe grounding in case of a wire fixture or strain relief failure. See drawings!*

Cable colors single-phase

Wiring code			
Color Code (IEC 60757)	Number Code	Use	Color
bn	2	L1 (phase)	brown
bu	1 (blue label)	N (neutral)	blue
Color Code (IEC 60757)	Bicolored		
gnye	GND (grounding)	green-yellow (protective grounding conductor)	

Note

Only replace defective cables with same type. Mind insulation type and class, wire cross-section, length, flexibility and color or number code!

Cable colors three-phase

Wiring code

Color Code (IEC 60757)	Number Code	Use	Color
bk	4	phase	black
bn	3	phase	brown
gy	2	phase	grey
bu	1 (blue label)	N (neutral)	blue

Color Code (IEC 60757)

Bicolored


gnye

GND (grounding)

green-yellow
(protective grounding conductor)

Note

Only replace defective cables with same type. Mind insulation type and class, wire cross-section, length, flexibility and color or number code!

Steps to finish the installation

Perform the following steps to finish the installation:

- Run an electrical safety test according to the local regulations if applicable. Mind the device's type plate for insulation class!
- Check external power switches and other installed safety devices for functionality.
- Inform the operator/customer about external power switches or the location of the power plug.
- Before switching on the device for the first time, see operation manual for information about how to establish suitable startup conditions!

 	Notice
	Risk of device damage by faulty operation • Consult operation manual for information about working liquids or other initial device preparations for a functional test

Customer service

Periodic testing of an installation

In many countries, industrial and commercial-building installations must be re-tested periodically by authorized agents. Please mind local laws and regulations!

Only authorized service personnel are allowed to perform repair work on the device. Authorization requires a comprehensive technical training and knowledge of possible dangers which might arise when working at the device. Such training and knowledge can only be provided by Buchi.

Addresses of official Buchi customer service offices are given on the Buchi website under: www.buchi.com. If malfunctions occur on your device or you have technical questions or application problems, contact one of these offices.

The customer service offers the following:

- Spare part delivery
- Repairs
- Technical advice

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BUCHI (Thailand) Ltd.
BÜCHI Labortechnik AG