



Application Note – N°. 904/2025

# On-Site sugar content analysis in cereal

**Abstract:** Breakfast cereals are a common source of nutrients but often contain high sugar levels, raising health concerns. Accurate sucrose measurement is vital for labeling and consumer awareness. This study assesses the ProxiScout™ handheld NIR scanner as a rapid, non-destructive, and cost-effective alternative to traditional HPLC. Results show ProxiScout™ delivers reliable sucrose estimates, supporting its use in quality control.



## 1. Introduction

---

Nutrition labels may not always reflect the actual sugar content of commercial products, with studies showing discrepancies of up to 20%. Given the health implications of excessive sugar consumption, there is a need for real-time, accurate methods to assess sugar levels in processed foods.

Traditional sugar analysis techniques, such as HPLC and enzymatic methods, are accurate but expensive and labor-intensive. Near-Infrared (NIR) spectroscopy presents a rapid, non-destructive alternative. This study demonstrates the effectiveness of ProxiScout™ in measuring sucrose levels in breakfast cereals, allowing for efficient on-site testing.

## 2. Equipment

---

The ProxiScout™ handheld NIR scanner was used to collect spectral data from breakfast cereal samples. The device enables real-time quality control and screening with minimal training, utilizing diffuse reflectance spectroscopy to capture spectral fingerprints.

- Setup: Scanner on Diffuse reflection mode coupled with the rotator accessory rotating stage for spectral averaging of Heterogeneous samples. (Figure 1).
- Spectral Range: 1350–2550 nm.



*Figure 1: ProxiScout™ Scanner coupled with the rotating accessory.*

## 3. Samples and Measurement Conditions

---

Samples:

- Total Samples: 60 breakfast cereal samples coated with sucrose, manufactured by a leading snack company in Ohio, USA.
- Sample Preparation:
  - Cereals were blended to achieve a homogeneous particle size.
  - Samples were transferred to glass petri dishes and placed on a rotating stage to ensure spectral averaging for reproducible measurements.

Measurement Conditions:

- Scan Time: 5 seconds per sample.
- Resolution: 16 nm at  $\lambda = 1,550$  nm.
- Spot Size: 10 mm<sup>2</sup>.
- Temperature: Room temperature.
- Averaging: Each sample was measured three times, with the results averaged for analysis.



## 4. Procedure

### 4.1 Reference Methods

- HPLC was used to determine sucrose levels in the samples.
- The sugars were separated using a stainless steel 7.8 mm ID x 300 mm Aminex® HPX-87C carbohydrate column under isocratic conditions at 80 °C using HPLC grade water with a flow rate of 0.6mL/min for 30 minutes.

### 4.2 Calibration Model Development

Partial Least Squares (PLS) regression models were developed to establish a linear relationship between spectral data and laboratory-measured sucrose content. The models were optimized using:

- Cross-validation techniques to enhance predictive accuracy.
- Generalized sample calibration, ensuring robustness against variations in sample composition.

### 4.3 Data Analysis

The PLS model was evaluated using cross-validation, which involved:

- Splitting data into calibration and validation sets.
- Train the model using the calibration set.
- Testing the model on the validation set to assess predictive performance.

The statistical parameters used for evaluation include:

- $R^2$  (Coefficient of Determination): Measures correlation between predicted and actual values (closer to 1 is better).
- RMSE (Root Mean Square Error): Represents prediction error (lower is better).
- SECV (Standard Error of Cross-Validation): Indicates model reliability (lower values suggest better accuracy).

## 5. Result

### 5.1 Model Performance Metrics

Sucrose Content Prediction:

- $R^2 = 0.99$ , indicating excellent correlation with HPLC results.
- SECV = 1.2%, demonstrating high accuracy.
- Performance surpasses previous NIR-based methods and bench-top FT-NIR systems.




Figure 2: Spectrum obtained from cereal samples taken with ProxiScout™-Module.




Figure 3: Partial Least Squares Regression (PLSR) plots for sucrose content in breakfast cereals by ProxiScout™-Module.

## 5.2 Key Observations

- The ProxiScout™ scanner accurately predicted sucrose concentrations in cereals.
- Minimal error compared to traditional HPLC analysis.
- The system enables real-time testing, reducing reliance on costly lab-based methods.
- The sucrose levels in the cereal samples ranged from 0 to 30 g/100 g, covering a wide range of concentrations resembling commercial breakfast cereals with reported averages of 22.5±12.6 g/100 g (Pombo-Rodrigues et al., 2017).
- PLSR Regression Results: SECV for sucrose was 1.2%, with a correlation coefficient (r) of 0.99.
- Comparison with other NIR Methods: ProxiScout™-Module demonstrated superior performance compared to bench-top FT-NIR (SEP = 1.5%) and another handheld analyzer (1.5%).
- NIR Spectra Findings: Characteristic absorption bands at 4390 cm<sup>-1</sup> (C-H bonds), 4420 cm<sup>-1</sup> (C-H in carbohydrates), 4880 cm<sup>-1</sup> (amide absorption in protein), 5180 cm<sup>-1</sup> (O-H in water), and 5807 cm<sup>-1</sup> (aliphatic C-H).
- Unique Sucrose Band: Found at ~4800 cm<sup>-1</sup>, confirming the presence of crystalline sucrose.

## 6. Conclusion

---

### 6.1 Comparison to Standard Methods

- ProxiScout™ results closely matched HPLC measurements, validating its suitability for on-site testing.
- Unlike laboratory methods, ProxiScout™ provides instant results, eliminating sample transfer and preparation delay.

The ProxiScout™ handheld NIR scanner has demonstrated excellent performance in measuring sucrose content in breakfast cereals. Its ability to provide rapid, non-destructive analysis makes it an ideal tool for cereal manufacturers, nutritionists, and regulatory agencies. By integrating real-time testing into production processes, ProxiScout™ enhances quality control and ensures compliance with nutrition labeling requirements.

## 7. References

---

- [2] Walker, R. W., & Goran, M. I. (2015). Labeling accuracy of added sugars in breakfast cereals and nutrition bars marketed to children. *British Journal of Nutrition*, 114(8), 1287-1294. <https://doi.org/10.1017/S0007114515002822>
- [3] Te Morenga, L., Mallard, S., & Mann, J. (2012). Dietary sugars and body weight: Systematic review and meta-analyses of randomized controlled trials and cohort studies. *BMJ*, 346, e7492. <https://doi.org/10.1136/bmj.e7492>