



## Recirculating Chiller F-100 / F-105

# Pre-installation Checklist

The following checklist will assist you in site preparation for your BUCHI instrument. Careful attention to the requirements presented in this document will help simplify instrument installation and ensure that the instrument operates as it was designed.

# Recirculating Chiller F-100 / F-105

## Pre-installation checklist

All appropriate items in this checklist must be completed, signed and returned to BUCHI. Return your checklist as soon as your laboratory is ready and functional for prompt installation of your recirculating chiller. If you are performing self-installation, please check the box for self-installation below and return the form to BUCHI.

### Environmental requirements

Ambient Temperature: 5 °C to 35 °C.

Humidity: Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C.

Altitude: max. 2000 m above sea level.

### Space requirements

Enough space as defined below: (Width x Depth x Height).

Min. 280 x 500 x 420 mm.

The instrument must have a clearance of 40 cm between itself and the wall, both in front and behind (sufficient cooling).

To ensure that the power can be cut by unplugging in case of an emergency, the mains plug must not be blocked by the instruments or any other items!

### Location of the unit

The Recirculating Chiller F-100 and F-105 must be installed on a clean, level, vibration free surface.

### Power requirements

The F-100 / F-105 requires its own standard mains outlet socket. A standard mains cable is supplied.

### Detailed Technical requirements

I have reviewed the detailed Technical Data on [www.buchi.com/en/support/downloads](http://www.buchi.com/en/support/downloads)  
My facility fulfills the required specifications.

# Recirculating Chiller F-100 / F-105

## Pre-installation checklist

### Facility sign-off

Return your checklist as soon as your facility is ready and functional for prompt installation of your Recirculating Chiller F-100 / F-105.

Environmental Requirements as specified for your facility.

Space Requirements.

Power Requirements.

Other Requirements.

Self-Installation.

### Please enter the following information and sign below:

Company name

Customer name

Email address

Purchase order n°.

I certify that the preparations required for the installation of the Recirculating Chiller F-100 / F-105 have been completed and functional as described in this document.

**I further acknowledge that if any of these requirements are not ready by the agreed upon installation date and I have not rescheduled the installation, then I will be responsible for any additional costs relating to the installation of the system.**

Customer signature

Date

Please sign and return by email to BUCHI as soon as your laboratory is ready for installation:

**Attention:** Find your BUCHI affiliate or partner on: [www.buchi.com/worldwide](http://www.buchi.com/worldwide)

**If you have any questions or queries, then please do not hesitate to contact us.**

# Core messages to our customers

## BUCHI creates added value

“Quality in your hands” is the guiding principle that shapes our philosophy and our actions. It challenges us to provide outstanding services that are precisely tailored to your needs. This means that we must stay in close contact with our customers. That is why we keep in touch and continue to work very hard to understand you and your business even better.

We help you by providing high-quality products, systems, solutions, applications and services that offer you added value. This allows you to focus entirely on your processes and your work.



### Competent

We have the technological expertise and decades of experience needed to provide competent support and work with you to continually improve our services.



### Reliable

We guarantee the quality and functionality of our equipment and will continue to help you quickly and efficiently whenever something does not operate to your satisfaction.



### Safe

By collaborating closely with you, we do everything in our power to make our products, systems, solutions, applications and services as safe as possible for people and the environment.



### Cost-effective

We strive to create a high level of economic benefit and maximum added value for you.



### Global

As an international family-owned business with own subsidiaries and qualified distributors, we have a presence wherever you are located.



### Easy

We support you by providing carefully designed solutions as well as instruments and systems that are easy to operate.



### Sustainable

We support environmentally friendly processes and manufacture products that have a long service life. We utilize advanced technologies to leave the smallest environmental footprint possible.

---

We are represented by more than 100 distribution partners worldwide.  
Find your local representative at:

[www.buchi.com](http://www.buchi.com)

Quality in your hands

