




# Fast, Accurate PFAS Recovery

## SyncorePlus® Polyvap

Improve PFAS sample concentration throughput, streamline your workflow, and provide accurate, reproducible reporting with the SyncorePlus® Polyvap. Benefit from:



**High Throughput**  
Concentrate up to 48 PFAS samples to final volume in parallel.

**Free Up Hood Space**  
Condensers enable high recovery of solvents to maintain a safe lab environment even outside the hood.

**Contamination-Free**  
Transport racks offer direct transfer from SPE to SyncorePlus® for workflow synergy. PFAS-free design with individual tube positions and vacuum ports prevent cross-contamination.



**Fast & Accurate Recovery**  
Automated, vacuum assisted evaporation ~2x faster than nitrogen blowdown. Low boiling point concentration more effectively preserves light-end PFAS.



**Improve Workflow**  
Use Interface Controller to automate process and enable unattended operation.

**Convenient Analytical Transfer**  
Reach final volume directly in 15 or 50-mL Falcon Tubes for easy transfer to HPLC, LC/MS, or combustion IC.


### Learn More:



On-Demand Webinar  
SyncorePlus® for PFAS



Stream  
SyncorePlus® Demo



Download  
Product Information



Ask an Expert  
Contact Us