

Certificate of Analysis

Kjeldahl Tablets Titanium

Order No.: 11072627
Product: Kjeldahl Tablets Titanium

Version created by: GOHH
Version creation date: 04.11.2023
Manufacturing date: 09.11.2023
Expiry date: 09.11.2028

Lot Number 04091123



Parameter	Specification	Result
Appearance	Round white with green/blue speckled flat tablets	Complies
Potassium sulfate	3.50 g/tablet	3.50 g/tablet
Copper sulfate pentahydrate	0.105 g/tablet	0.105 g/tablet
Titanium dioxide	0.105 g/tablet	0.105 g/tablet
Weight uniformity	3.71 ± 0.2 g	3.70 ± 0.02 g
Friability	≥ 100 N/mm ²	244 ± 50 N/mm ²
Nitrogen	≤ 12 ppm	≤ 6 ppm

Storage: Tightly closed in dry and well-ventilated rooms

This document has been generated electronically and is valid without a signature.

BÜCHI Labortechnik AG

Meierseggstrasse 40
9230 Flawil
Switzerland

T +41 71 394 63 63
F +41 71 394 64 64

www.buchi.com

Quality in your hands